Assay for the Qualitative Detection of Alcohol in Urine or Saliva

Part No. 47464
English Intended Use

On-Site Alcohol® is intended for the rapid detection of alcohol in saliva or urine. On-Site Alcohol is recommended for professional use only and is not intended for over-the-counter sales to the general public. Professional use includes applications in health care, corrections, and drug treatment.

Summary and Explanation of Test


Alcohol is the most commonly abused drug, and 200,000 Americans die every year of alcohol abuse. Five to ten percent of the US population is dependent on alcohol; alcohol abuse is the leading killer of Americans aged 15-45. Alcohol abuse is also a frequent contributing factor in all types of accidents and criminal violations.

In several studies, a significant correlation between saliva and blood ethanol concentrations has been found.¹ The saliva/blood ratio is approximately 1.07 when both fluids are collected simultaneously.¹ However, due to metabolic and physiologic variables, this correlation is not always constant. The analysis of urine as an indirect means of determining an equivalent blood alcohol concentration is a long established technique, but the relationship between the concentrations of alcohol in urine and blood may vary widely.² Confirmation of On-Site Alcohol results using gas chromatography on whole blood is recommended if the specific blood alcohol concentration level must be determined.

Test Principle

The On-Site Alcohol assay is based on two coupled enzymatic reactions involving alcohol dehydrogenase (ADH), the coenzyme nicotinamide adenine dinucleotide (NAD), and diaphorase.

In the first reaction step, ethanol in the sample is oxidized to acetaldehyde by ADH. The coenzyme NAD is concomitantly reduced to NADH. The second step uses the enzyme diaphorase to catalyze the oxidation of NADH to NAD. The hydrogen resulting from this reaction is transferred to a colorless tetrazolium salt to produce a highly colored formazan dye.³

In practice, urine or saliva specimens are deposited on a chemically treated pad that volatilizes reducing alcohols. The resulting vapors are concentrated onto a membrane saturated with enzyme solution and tetrazolium salt. When the enzymetetrazolium salt solution is exposed to ethanol vapor, a formazan dye is formed, producing a purple positive sign (+) in the center of the pale yellow membrane.

Reagents

Each On-Site Alcohol kit contains 50 individually packaged alcohol test cards, one alcohol test per package.

Each package contains:

1. Test Card: One plastic test card with reaction site for testing one specimen.
2. Alcohol Detection Reagent: 75 µL of tetrazolium dye solution contained in the detection reagent well.
3. Disposable Transfer Pipet (Small): Small plastic pipet used to transfer Detection Reagent from the detection reagent well to the detection reagent pad on the test card.
4. Disposable Swab: Absorbent swab used to transfer saliva specimen directly to the test card.

5. Disposable Transfer Pipet (Large): Large plastic pipet used to transfer urine or saliva specimen to the test card.
6. Desiccant: One packet containing a drying substance to ensure the test card remains dry.

Precautions and Warnings

Exercise the normal precautions for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines.

1. The test card is moisture-sensitive and must be used immediately after opening the package.
2. Avoid cross-contamination by using new pipets for each specimen.

Storage and Stability

1. Test cards can be stored at 59 - 86 °F (15 - 30 °C) with factory seals intact.
2. Do not allow test reagent to freeze.
3. Do not use kit beyond the expiration date that is printed on individual On-Site Alcohol packages. If foil is damaged (hole or tear), do not use kit.

Specimen Collection and Handling

Collect specimens in clean, sealable containers. Keep the containers closed tightly when not in use. Test the specimens as soon as possible after collection.

Specimens containing human-sourced materials should be handled as if potentially infectious. Use safe laboratory procedures such as those outlined in *Biosafety in Microbiological and Biomedical Laboratories* (HHS Publication Number [CDC] 93-8395).

Saliva may be collected in a container or directly onto the disposable swab without the use of a specimen collection container. (See Step 3, Method A) If a container is used, be sure to minimize the air volume above the specimen. If a swab is used, the swab must be immediately inserted into the prepared On-Site Alcohol test card. Note: After the test subject eats or drinks, it is important to wait at least 10 minutes before collecting saliva.

Procedure
Materials Provided

50 individually packaged alcohol test cards, one test per package.

Ancillary Materials

	Part No.	Quantity
Gloves	42210	50 pairs
Collection cups	47123	100 cups/box
Adulteration strips	3111407	25 strips/box
Evidence tape	42212	100 strips/pkg

Instructions for Use

Note: The On-Site Alcohol test should be performed at room temperature, 65 - 85 °F (18 - 29 °C).

Step 1: Open the On-Site Alcohol test card by peeling the lid from left to right. Discard the desiccant packet. Record specimen identification in the designated area marked "I.D." on the front of the test card.

Step 2: Detection Reagent Preparation

Use the small disposable transfer pipet provided to transfer one drop of detection reagent from the detection reagent well to the detection reagent pad located to the right of the well. To do so, squeeze the bulb end of the pipet, place the tip in the detection reagent well, and release the bulb to draw the detection reagent into the pipet. Hold the pipet vertically over the detection

reagent pad to slowly dispense one drop. Use one drop only. Holding the pipet vertically (90° angle to the card) will help reduce the occurrence of air bubbles.

WARNING: Avoid Detection Reagent contact with eyes, skin and mucous membranes. If contact occurs, rinse immediately and thoroughly with water. **KEEP OUT OF REACH OF CHILDREN.**

Step 3: Specimen Application-Use Method A or B.

Method A

Saliva Testing: Use the disposable swab for direct saliva specimen collection following Detection Reagent preparation (Step 2). Instruct the test subject to completely saturate the disposable swab with saliva by pressing the foam head against the inner mouth area and tongue to absorb as much saliva as possible in 10 to 20 seconds. Transfer the saturated swab to the test card and insert the swab vertically into the Specimen well. The swab will press-fit into the specimen well. The swab must remain in this position for the duration of the test. Note: For reliable results, the disposable swab foam head must be thoroughly saturated with saliva. If saliva is collected in a tube, it can be transferred to the Specimen well via the large pipet provided. Note: Due to its viscosity, one discrete drop of saliva may be difficult to dispense from the pipet. In this case, fill the Specimen well approximately one quarter full of saliva from the pipet.

Method B

Urine Testing: Use the large disposable pipet provided to transfer one drop (45 µL) of urine from the specimen container into the Specimen well. To do so, squeeze the bulb end of the pipet, place the tip in the specimen, and release the bulb to draw up a small amount of specimen into the pipet. Hold the pipet vertically over the Specimen well, allowing the drop to form and fall into the well. Gently squeeze the bulb to dispense one drop only. Holding the pipet vertically (90° angle to the well) will help reduce the occurrence of air bubbles.

Step 4: Lay the test card flat and wait two minutes before reading and recording the results. Read the test results at two minutes to ensure On-Site Alcohol testing accuracy. Note: Dispose of used card and components properly. Keep out of reach of children.

Quality Control

In determining and maintaining assay reliability and performance, the use of both positive and negative controls at the start of each test or group of tests is recommended. Distilled or deionized water can be used for a negative control.

These controls should be run to ensure proper technique and assay reliability and performance. If the expected results are not achieved with the positive and negative controls, do not proceed with testing. Contact the Varian Technical Support at 1-800-737-9667, for technical assistance or information about commercially-available controls.

Results

Negative - Detection reagent pad remains pale yellow for two minutes.

Preliminary Positive - A purple "positive" sign (+) is produced in the center of the pale yellow detection pad in two minutes or less. A purple + sign, independent of intensity, produced in the center of the detection reagent pad, is recorded as a presumptive positive. An incomplete or broken + sign should be interpreted as a presumptive positive. On-Site Alcohol results with saliva provide a means of estimating whether blood alcohol concentration is greater than 0.02% w/v (20 mg/dL weight/volume ratio). Presumptive positive On-Site results with urine are indicative of recent alcohol ingestion but do not necessarily indicate intoxication, impairment or blood alcohol concentrations over 0.02% w/v (20 mg/dL). Confirmation of On-Site Alcohol results using gas chromatography on whole blood is recommended if blood alcohol concentration must be determined.

Procedural Limitations

On-Site Alcohol is a qualitative test for the detection of ethanol. Methanol and isopropanol will not produce positive results if present in urine or saliva. Patients with presumptive positive results who demonstrate euphoria or in whom toxic levels are suspected should be further tested

using a quantitative whole blood method, preferably gas chromatography. Capillary and arterial blood samples most accurately reflect brain ethanol concentrations.

Care should be taken to avoid contamination of specimens and equipment. The atmosphere must be free of alcohol or other solvent vapors, as they may affect results. Areas of high humidity may shorten shelf life.

Expected Values

Alcohol is not normally present in human urine or saliva in measurable quantities. The total amount of ingested alcohol eliminated in the urine varies depending on the rate of absorption and the volume of urine excreted. In the post-absorptive stage, urine-alcohol concentration will be approximately 25% higher than blood-alcohol levels.⁴ Saliva (salivary ethanol) has been shown to be a good index of blood-ethanol concentration provided the saliva is obtained at least 20 minutes after the ingestion of alcohol.⁵ Persons with severe ethanol intoxication can have urine and blood ethanol concentrations ranging as high as 0.8% w/v (800 mg/dL) ethanol.

Performance Characteristics

Accuracy

In clinical trials, On-Site Alcohol results were compared to results obtained using the (previous) acid-dichromate-based On-Site Alcohol assay, which has shown accuracy of 98% in, correlation studies with gas chromatography and accuracy of 100% when correlated with a quantitative enzymatic saliva method.

Using data from controlled (internal) and field site (external) testing (see Summary of Data Table), a correlation of 97% (180/186) was obtained for the acid-dichromate-based On-Site Alcohol assay and the new, enzymatic version. In the controlled tests, both the enzymatic and the acid-dichromate-based versions were used to test spiked samples (16), specimens from non-drinking volunteer subjects (20), and specimens from volunteers who had consumed measured quantities of ethanol (58). The drinking volunteer subjects donated saliva and urine 30 minutes to one hour after the alcohol consumption. In some cases, insufficient alcohol was consumed to produce concentrations 0.02% w/v, and these test results were defined as negative.

Summary of Data

	Positives Tested		Positive Results		Negatives Tested		Negative Results	
	Urine	Saliva	Urine	Saliva	Urine	Saliva	Urine	Saliva
Internal	34	34	35	34	11	14	10	14
External*	7	3	9	4	25	58	23	57
Subtotal	41	37	44	38	36	72	33	71
Totals	78		82		108		104	

Positives are defined as samples that tested positive by On-Site Alcohol at $\geq 0.02\%$ w/v ethanol.

Negatives are defined as samples that tested negative by On-Site Alcohol at $\leq 0.02\%$ w/v ethanol.

*External studies were conducted at four different probation/parole testing facilities.

Accuracy was further demonstrated by running fifty-two (52) replicates of three ethanol controls. After application of the sample, the test was timed and after exactly two minutes, a reading was taken. Results were reported as either positive or negative.

	Positives	Negatives	% Reported as Positive
Control 1 (0.008% w/v)	0	52	0
Control 2 (0.020% w/v)	36	16	69
Control 3 (0.032% w/v)	50	2	96

Precision

Precision was determined by running fifty-two (52) replicates of three ethanol controls. The time necessary to produce an observable positive following the application of the control sample was carefully observed and recorded for each test. The mean, standard deviation, and coefficient of variation were calculated based on the timed test reactions.

	Mean	Standard Deviation	Coefficient of Variation
Control 1 (0.008% w/v)	155 sec (2 min 35 s)	15 sec	9.5%
Control 2 (0.020% w/v)	118 sec (1 min 58 s)	9 sec	8.0%
Control 3 (0.032% w/v)	108 sec (1 min 48 s)	10 sec	9.0%

Sensitivity

Twenty replicates of six ethanol controls were run to determine sensitivity. After application of the samples, the test was timed and after exactly two minutes a reading was taken. Results were reported as either positive or negative.

	Positives	Negatives	% Reported as Positive
Control 1 (0.005% w/v)	0	20	0
Control 2 (0.010% w/v)	0	20	0
Control 3 (0.015% w/v)	11	9	55
Control 4 (0.020% w/v)	19	1	95
Control 5 (0.025% w/v)	20	0	100
Control 6 (0.030% w/v)	20	0	100

Specificity

On-Site Alcohol is designed to detect ethyl alcohol (ethanol). The following substances were tested with On-Site Alcohol and were found to be non-interfering at the stated concentrations.⁶ The only exception is 1-propanol, which produced a faint positive with On-Site Alcohol at a concentration of 0.05% w/v. It is unlikely that such a concentration of 1-propanol would be found in human urine or saliva.

Substance	Concentration (% w/v)	Result
Acetaldehyde	1.00%	neg
Acetone	1.00%	neg
1-Butanol	0.10%	neg
2-Butanol	1.00%	neg
Ethylene glycol	1.00%	neg
Glycerol	1.00%	neg
Isopropanol	0.10%	neg
Methanol	1.00%	neg
1-Propanol	0.05%	pos

Caution: Ingestion or general use of over-the-counter medications and products containing alcohol can produce positive specimens when tested with On-Site Alcohol. Vicks[®], Nyquil[®] Nighttime Cold Medicine, Binaca[®] Concentrated Breath Spray, and Listerine[®] Antiseptic are examples of these products. Subjects claiming to have recently used such products should be observed continuously for 20 minutes and retested.

NOTICE:

On-Site Alcohol[®] Testing for National Highway Traffic Safety Administration (NHTSA), Department of Transportation (DOT), Highway and Safety Programs:

On-Site Alcohol has been designed for use at room temperature (65 - 85 °F or 18 - 29 °C). This includes uses in health care, corrections, and drug treatment, where tests are generally performed indoors at controlled temperatures. In screening for alcohol in bodily fluids under NHTSA DOT Highway Safety Programs, it is anticipated that On-Site Alcohol tests may be performed outdoors at varying temperatures. It has been determined that On-Site Alcohol will maintain its performance characteristics under a range of environmental conditions when adjustment is made in the time the result is read (see below). Testing at temperatures outside this range will not produce acceptable results.

Temperature Range

Low	Normal	High
50 - 64 °F (10 - 17 °C)	65 - 85 °F (18 - 29 °C)	86 - 104 °F (30 - 40 °C)

Read Result at:

five minutes after specimen application	two minutes after specimen application	one minute after specimen application
---	--	---

References

1. Jones, A.: Distribution of Ethanol Between Saliva and Blood in Man. *Clin. and Exper. Pharmacol. & Physiol.* 6, 53-59, 1979.
2. Kaye, S., and Cardona, E.: Errors of Converting a Urine Alcohol Value into a Blood Alcohol Level. *Amer. J. of Clin. Path.*, 52(5), 577-584, 1969.
3. Bergmeyer, H.: *Methods of Enzymatic Analysis*, Academic Press: New York, p. 285, 1965.
4. *Review of Medical Pharmacology*, 4th ed. ; Meyers, P. , Jawetz, E., and Goldfien, A., Eds.; Lange Medical: Los Altos, CA, pp 236-238, 1974.
5. McColl, K., Whiting, B., Moore, M., and Goldberg, A.: Correlation of Ethanol Concentration in Blood and Saliva. *Clinical Science*, 56, 283-286, 1979.
6. Data on file, Varian, Inc.

Urine or Saliva (using the large transfer pipet) Testing Procedure

1 Using the small transfer pipet, pick up reagent from the detection reagent well at the left side of the test card.

2 Transfer the reagent from the detection reagent well to the detection reagent pad (Result). Add only one drop of reagent to the well.

3 Use the large disposable pipet to pick up the urine or saliva specimen from the specimen container.

4 Transfer one drop of the urine or saliva specimen to the specimen well (Sample) of the test card.

5 Read results 2 minutes after specimen application. A purple positive sign in the detection reagent pad (Result) at ≤ 2 minutes indicates ethanol concentration of $\geq 0.02\%$ w/v. Negative specimens ($<0.02\%$) do **not** produce a "+" sign in ≤ 2 minutes.

Saliva Testing Procedure

1 Using the small transfer pipet, pick up reagent from the detection reagent well at the left side of the test card.

2 Transfer the reagent from the detection reagent well to the detection reagent pad (Result). Add only one drop of the reagent to the well.

3 Use the disposable swab for direct saliva collection. Instruct the test subject to press the foam head against the inner mouth area and tongue for 10 to 20 seconds, or until the swab is saturated.

4 Transfer the saturated swab to the test card, and insert the swab vertically into the specimen well (Sample). The swab must remain in this position for the duration of the test.

5 Read results 2 minutes after specimen application. A purple positive sign in the detection reagent pad (Result) at ≤ 2 minutes indicates ethanol concentration of $\geq 0.02\%$ w/v. Negative specimens ($< 0.02\%$) do **not** produce a "+" sign in ≤ 2 minutes.

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Literatur

1. Jones, A.: Distribution of Ethanol Between Saliva and Blood in Man. *Clin. and Exper. Pharmacol. & Physiol.* 6, 53-59, 1979.
2. Kaye, S., and Cardona, E.: Errors of Converting a Urine Alcohol Value into a Blood Alcohol Level. *Amer. J. of Clin. Path.*, 52(5), 577-584, 1969.
3. Bergmeyer, H.: *Methods of Enzymatic Analysis*, Academic Press: New York, p. 285, 1965.
4. *Review of Medical Pharmacology*, 4th ed.; Meyers, P., Jawetz, E., and Goldfien, A., Eds.; Lange Medical: Los Altos, CA, pp 236-238, 1974.
5. McColl, K., Whiting, B., Moore, M., and Goldberg, A.: Correlation of Ethanol Concentration in Blood and Saliva. *Clinical Science*, 56, 283-286, 1979.
6. Dokumentation, Varian, Inc.

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Bibliographie

1. Jones, A.: Distribution of Ethanol Between Saliva and Blood in Man. Clin. and Exper. Pharmacol. & Physiol. 6, 53-59, 1979.
2. Kaye, S., and Cardona, E.: Errors of Converting a Urine Alcohol Value into a Blood Alcohol Level. Amer. J. of Clin. Path., 52(5), 577-584, 1969.
3. Bergmeyer, H.: *Methods of Enzymatic Analysis*, Academic Press: New York, p. 285, 1965.
4. *Review of Medical Pharmacology*, 4th ed. ; Meyers, P. , Jawetz, E., and Goldfien, A., Eds.; Lange Medical: Los Altos, CA, pp 236-238, 1974.
5. McColl, K., Whiting, B., Moore, M., and Goldberg, A.: Correlation of Ethanol Concentration in Blood and Saliva. *Clinical Science*, 56, 283-286, 1979.
6. Documentation de Varian, Inc.

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Riferimenti

1. Jones, A.: Distribution of Ethanol Between Saliva and Blood in Man. *Clin. and Exper. Pharmacol. & Physiol.* 6, 53-59, 1979.
2. Kaye, S., and Cardona, E.: Errors of Converting a Urine Alcohol Value into a Blood Alcohol Level. *Amer. J. of Clin. Path.*, 52(5), 577-584, 1969.
3. Bergmeyer, H.: *Methods of Enzymatic Analysis*, Academic Press: New York, p. 285, 1965.
4. *Review of Medical Pharmacology*, 4th ed. ; Meyers, P. , Jawetz, E., and Goldfien, A., Eds.; Lange Medical: Los Altos, CA, pp 236-238, 1974.
5. McColl, K., Whiting, B., Moore, M., and Goldberg, A.: Correlation of Ethanol Concentration in Blood and Saliva. *Clinical Science*, 56, 283-286, 1979.
6. Dati presenti nel file, Varian, Inc.

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Bibliografia

1. Jones, A.: Distribution of Ethanol Between Saliva and Blood in Man. *Clin. and Exper. Pharmacol. & Physiol.* 6, 53-59, 1979.
2. Kaye, S., and Cardona, E.: Errors of Converting a Urine Alcohol Value into a Blood Alcohol Level. *Amer. J. of Clin. Path.*, 52(5), 577-584, 1969.
3. Bergmeyer, H.: *Methods of Enzymatic Analysis*, Academy Press: New York, p. 285, 1965.
4. *Review of Medical Pharmacology*, 4th ed. ; Meyers, P., Jawetz, E., and Goldfien, A., Eds.; Lange Medical: Los Altos, CA, pp 236-238, 1974.
5. McColl, K., Whiting, B., Moore, M., and Goldberg, A.: Correlation of Ethanol Concentration in Blood and Saliva. *Clinical Science*, 56, 283-286, 1979.
6. Ficheiros da Varian, Inc.

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© 2003, Varian, Inc.
Varian B.V.