English

Intended Use

OnTrak TesTstik™ for Amphetamines is an *in vitro* diagnostic test intended for professional use for the qualitative detection of amphetamines in urine at or above a cutoff concentration of 1000 ng/mL. OnTrak TesTstik for Amphetamines is not intended for over-the-counter sale.

OnTrak TesTstik provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmation method.¹ Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Summary

The amphetamines (amphetamine, methamphetamine and the structurally similar “designer” drugs, e.g., “Ecstasy”), are central nervous system (CNS) stimulants related chemically and pharmacologically to naturally occurring epinephrine and norepinephrine, the body’s catecholamines. The amphetamines increase wakefulness and suppress appetite through CNS stimulation, at the same time causing a sense of increased energy, self-confidence, well-being and euphoria.^{2,3} Cardiovascular effects include increased blood pressure and cardiac arrhythmias.³ Overdoses can result in irrational behavior, reduced mental acuity and psychosis.

Abuse patterns usually begin with oral drug use but may progress to intravenous administration culminating in “speed runs” of repetitive, sustained use with accompanying wakefulness, tension, and paranoia interspersed with “crashes” of deep sleep.⁴ The amphetamines are readily absorbed from the stomach following ingestion and rapidly enter the blood.⁵ As CNS stimulants, their effects include restlessness, anorexia, insomnia, agitation, tremors, and increased motor activity.^{2,3,5} They are metabolized by the liver by ring hydroxylation, demethylation (methamphetamine is demethylated to amphetamine), deamination, and conjugation.^{6,7} Substantial fractions of the amphetamines are excreted unchanged in the urine along with their metabolites; urinary pH affects the percentage of unchanged drug excreted.^{6,8}

Test Principle

OnTrak TesTstik assays are based on the principles of microparticle capture inhibition. The test relies on the competition between drug, which may be present in the urine being tested, and drug conjugate immobilized on membrane for binding to antibody coated onto colored microparticles. When the TesTstik is immersed in the urine specimen, some of the specimen is absorbed into the TesTstik sample pad. The absorbed specimen travels through a reagent strip contained in the device by capillary action. In the reagent strip, the specimen rehydrates and mobilizes antibody-coated blue microparticles. The microparticle-urine suspension continues to migrate through the reagent strip and comes in contact with the immobilized drug conjugate. In the absence of drug in the urine, the antibody-coated microparticles bind to the drug conjugate and a blue band is formed at the result window (“negative” sign).

When drug is present in the specimen, it binds to the antibody-coated microparticles. If sufficient drug is present, the microparticles are inhibited from binding the drug conjugate and no blue band is formed at the result window. A presumptive positive (“non-negative”) sample causes the membrane to remain white (“positive” sign).

An additional antibody/antigen reaction occurs at the “TEST VALID” area. The “TEST VALID” blue band forms when anti-bovine serum albumin (anti-BSA) antibodies, which are imbedded in the reagent

membrane, bind to the BSA antigen on the blue microparticles. The presence of the “TEST VALID” band indicates that the test has completed, the reagents are viable, and the results are ready to interpret.

Reagents

1. Blue-dyed microparticles coated with mouse monoclonal anti-amphetamine analogue and BSA in a buffered solution containing preservative and dried onto a membrane.
2. Drug conjugates immobilized on membrane.
3. Mouse monoclonal anti-BSA antibody immobilized on membrane.

Precautions and Warnings

For *in vitro* diagnostic use. Exercise the normal precautions for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines.

Storage and Stability

OnTrak TesTstik should be stored at room temperature, 15 - 30 °C (59 - 86 °F), in the original foil pouch. If the foil pouch containing the TesTstik is damaged (e.g., a hole or tear), do not use that TesTstik. Do not use the TesTstik after the expiration date printed on the foil pouch.

Specimen Collection and Preparation

OnTrak TesTstik is formulated for use with urine specimens. Fresh urine specimens do not require any special handling or pretreatment. No additives or preservatives are required. Specimens should be stored refrigerated (2 - 8 °C or 36 - 46 °F) and tested within two days of collection. For prolonged storage, freezing is recommended.

Specimens containing human-sourced materials should be handled as if potentially infectious. Use safe laboratory procedures such as those outlined in *Biosafety in Microbiological and Biomedical Laboratories* (HHS Publication Number [CDC] 93-8395).

Testing Procedure

Materials provided

50 individually packaged OnTrak TesTstik for Amphetamines sticks

Ancillary materials

	Part No.	Quantity
Evidence Tape	42212	100 strips/pkg.
Latex Gloves	42210	50 pairs
Collection Cups	47123	100 cups/box
Adulteration Strips	3111407	25 strips/box

Contact Customer Technical Support for information regarding commercially available control material.

Assay

OnTrak TesTstik is ready for use directly from the foil pouch and must be used within eight hours after removal from the pouch. This time frame may be shorter in high humidity climates. Do not open the foil pouch unless it is at room temperature (15 - 30 °C or 59 - 86 °F).

Specimens should be between room temperature and body temperature 15 - 37 °C (59 - 99 °F). If the specimen was refrigerated, allow it to reach room temperature before testing. To ensure a sufficient quantity of sample for confirmation testing (as necessary), collect a minimum of 30 mL of sample.

1. To expose the sample pad, pull the protective slide back until it stops (up to the “Dip Line”).
2. Gently immerse the exposed sample pad into the specimen until the specimen reaches the dip line. Do not immerse the entire device into the specimen. Keep the OnTrak TesTstik for Amphetamines in the specimen for approximately 10 seconds.

- Remove the OnTrak TesTstik for Amphetamines from the specimen and while holding the OnTrak TesTstik for Amphetamines over the specimen, push the protective slide back over the sample pad until it stops. Allow the test to proceed until a blue band is observed in the "TEST VALID" window (see Results section below). The OnTrak TesTstik for Amphetamines may be placed down or held while the test is developing. There is no required timing step. Note: If the Test Valid band does not appear, re-dip the OnTrak TesTstik for Amphetamines for approximately 10 seconds. Allow to sit undisturbed for an additional 1-3 minutes. Once the Test Valid band forms, proceed to step 4.
- Remove the tab that covers the result window from the protective slide. Read and record the result.

Quality Control

Quality control testing at regular intervals is good laboratory practice and may be required by federal, state, or local guidelines. Always check with the appropriate licensing or accrediting bodies to ensure your quality control program meets the established standards.

Internal control: Each OnTrak TesTstik AMP has built-in process controls. The "TEST VALID" band is an internal POSITIVE process control. A blue "TEST VALID" band should always appear if an adequate sample volume is used, the sample and reagent are wicking on the membrane, and the reagents in the "TEST VALID" band and the conjugate-color indicator are working. In addition, the background in the result window should become clear and provide a distinct result. If smears of color appear in the background, nonspecific binding of the antibody may be occurring, making the test result difficult to interpret. The clearing of the result window background may be considered an internal NEGATIVE process control. If the "TEST VALID" band does not appear in the test window, or the background does not clear, then the test is not valid and a new test must be performed.

External control: The use of external positive and negative controls is recommended to test each shipment of product or when a new lot is used. Contact Customer Technical Support for information regarding commercially available control material. The controls should produce the expected results. If expected results are not obtained, do not run test specimens, and contact Customer Technical Support. Follow the appropriate federal, state and local guidelines when running external controls.

Results

Once a full, partial, or broken blue band is formed in the "TEST VALID" window (in approximately 3 minutes) remove the cover tab and interpret results as either presumptive positive or negative. Results may be interpreted for up to 30 minutes after the "TEST VALID" band forms.

A negative result is the presence of a blue band ("negative" sign) anywhere in the test result window. Since a negative sample may give a faint or incomplete color band, **any** full, partial, or broken band in the test result window indicates a negative result.

A presumptive positive result (drug present at or above the cutoff) is the lack of a blue band in the test result window. The window appears as a white plus (+) sign.

Limitations

See Specific Performance Characteristics for information on substances tested for cross-reactivity in this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors). Adulteration of the sample with oxidizing agents and/or acidic or alkaline agents can cause erroneous results. If adulteration is suspected, another sample should be collected and a new OnTrak TesTstik used. For information on how adulterants may affect the performance of this assay, contact Customer Technical Support.

A presumptive positive result with this assay indicates the presence of drug or drug metabolites in urine but does not reflect the degree of intoxication.

Specific Performance Characteristics⁹

Accuracy

OnTrak TesTstik for Amphetamines was evaluated using specimens screened by an automated immunoassay and confirmed positive by GC/MS (at a 500 ng/mL cutoff). Fifty samples positive for amphetamines were positive by OnTrak TesTstik (100%).

		GC/MS	
		+	
TesTstik	+	50	
	-	0	

One hundred six urine samples, obtained from a clinical laboratory and screened negative by an automated immunoassay relative to a 1000 ng/mL cutoff for amphetamines, were evaluated using OnTrak TesTstik. One hundred five were negative for amphetamine by OnTrak TesTstik (>99%). The one positive sample did not contain amphetamine by GC/MS. All positive and negative samples were also assayed by, and compared to, Abuscreen OnTrak for Amphetamines (1000 ng/mL). One hundred fifty six samples tested by both OnTrak TesTstik and Abuscreen OnTrak for Amphetamines demonstrated 97.4% agreement.

		Abuscreen OnTrak	
		+	-
TesTstik	+	48	3*
	-	1**	104

* Two samples contained 895 and 1282 ng/mL amphetamine, respectively, by GC/MS. One sample did not contain amphetamine by GC/MS.

**Sample did not contain amphetamine by GC/MS.

Precision

The precision of OnTrak TesTstik for Amphetamines was determined by testing various concentrations of drug standards for the primary drugs or drug metabolites detected by this system. All lots are required to minimally perform with a greater than 95% confidence level that a negative result will be attained with drug at 25% ($\frac{1}{4}x$) of its cutoff concentration and a greater than 95% confidence level that a positive result will be attained with drug at 150% ($1\frac{1}{2}x$) of its cutoff concentration.

Specificity

The following structurally similar compounds for amphetamines were tested for cross-reactivity with OnTrak TesTstik for Amphetamines. Each compound tested was prepared in normal human urine. The results are expressed as that amount of the compound capable of giving a result equivalent to the cutoff.

Amphetamine-related compounds	Approx. ng/mL equivalent to 1000 ng/mL (1:1) d,l-amphetamine	Approx. percent cross-reactivity
d,l-Amphetamine	1000	100
p-Hydroxyamphetamine	2000	50
Methylenedioxyamphetamine (MDA)	4000	25
l-Amphetamine	12,500	8
Phenylpropanolamine	>500,000	<0.2

Cross-Reactivity with Additional Drugs

The following compounds were tested at 100,000 ng/mL and found not to cross-react with the OnTrak TesTstik for Amphetamines assay.

Acetaminophen	Codeine
Acetylsalicylic acid	Dextromethorphan
Aminopyrine	Dextropropoxyphene
Amitriptyline	Diazepam
Amobarbital	Diphenhydramine
Ampicillin	Diphenylhydantoin
Ascorbic acid	Dopamine
Aspartame	Doxylamine
Atropine	Ecgonine
Benzocaine	Ecgonine methyl ester
Benzoyllecgonine	Ephedrine
Benzphetamine	Epinephrine
Butabarbital	Erythromycin
Caffeine	Estriol
Calcium hypochlorite	Fenoprofen
Chlordiazepoxide	Furosemide
Chloroquine	Gentisic acid
Chlorpheniramine	Glutethimide
Chlorpromazine	Guaiaicol glycerol ether
Cocaine	Hydrochlorothiazide

1 Gently tear open the foil pouch to remove the stick. Pull the protective slide back to expose the sample pad.

2 Immerse the sample pad into the specimen until the specimen reaches the DIP LINE. Hold the stick in the specimen for approximately 10 seconds.

3 Remove the stick from the specimen. While holding the stick over the specimen cup, push the protective slide back over the sample pad.


Test complete—
ready to interpret results


Test NOT complete—
DO NOT interpret results

4 Allow the test to proceed until a blue band forms in the “TEST VALID” window (see example on the LEFT). Test is NOT complete until a blue band forms in the “TEST VALID” window (see example on the RIGHT). Once the “TEST VALID” line appears, remove the cover tab and interpret results.

5 A blue band in the result window indicates a NEGATIVE result.

6 No color in the result window indicates a POSITIVE result. The window appears as a white plus (+) sign.

3-Hydroxytyramine
Ibuprofen
Imipramine
Isoproterenol
Ketamine
LSD
Lidocaine
MDMA
Melanin
Meperidine
Methadone
Methamphetamine
Methaqualone
Methylphenidate
Methyprylon
Morphine Sulfate
Naloxone
Naltrexone
Naproxen
Niacinamide
Nordiazepam
Norethindrone
Norpseudoephedrine

Oxazepam
Penicillin G
Pentobarbital
Phencyclidine
 β -Phenethylamine
Phenobarbital
Phenothiazine
Phentermine
Phenylbutazone
Phenylephrine
Procaine
Promethazine
Pseudoephedrine
Quinidine
Quinine
Secobarbital
Sulindac
 Δ^9 -THC-9-carboxylic acid
Tetracycline
Tetrahydrozoline
Trifluoperazine
Tyramine
Verapamil

References

1. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
2. Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc: 1984.
3. Hoffman B, Lefkowitz R. Catecholamines and Sympathomimetic Drugs. In: Gilman A, Rall T, Neis A, Taylor P, eds. *Goodman and Gillman's The Pharmacological Basis of Therapeutics*. 8th ed. New York, NY: McGraw Hill Pub Co; 1993:187-220.
4. Rowland M. Amphetamine blood and urine levels in man. *J Pharm Sci*. 1969;58(4):508-509..
5. Dimijian GG. Contemporary drug abuse. In: Goth A. ed. *Medical Pharmacology: Principles and Concepts*. 11th ed. St. Louis, Mo: CV Mosby Company; 1984:336-365..
6. Dring LG, Smith RL, Williams RT. The metabolic fate of amphetamine in man and other species. *Biochem J*. 1970;116:425-435.
7. Caldwell J. The metabolism of amphetamines in mammals. *Drug Metab Rev*. 1976;5(2):219-280.
8. Beckett AH, Rowland M. Urinary excretion kinetics of amphetamine in man. *J Pharm Pharmacol*. 1965;17:628-639
9. Data on file at Varian, Inc.

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Literatur

1. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
2. Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc: 1984.
3. Hoffman B, Lefkowitz R. Catecholamines and Sympathomimetic Drugs. In: Gilman A, Rall T, Neis A, Taylor P, eds. *Goodman and Gillman's The Pharmacological Basis of Therapeutics*. 8th ed. New York, NY: McGraw Hill Pub Co; 1993:187-220.
4. Rowland M. Amphetamine blood and urine levels in man. *J Pharm Sci*. 1969;58(4):508-509..
5. Dimijian GG. Contemporary drug abuse. In: Goth A. ed. *Medical Pharmacology: Principles and Concepts*. 11th ed. St. Louis, Mo: CV Mosby Company; 1984:336-365.
6. Dring LG, Smith RL, Williams RT. The metabolic fate of amphetamine in man and other species. *Biochem J*. 1970;116:425-435.
7. Caldwell J. The metabolism of amphetamines in mammals. *Drug Metab Rev*. 1976;5(2):219-280.
8. Beckett AH, Rowland M. Urinary excretion kinetics of amphetamine in man. *J Pharm Pharmacol*. 1965;17:628-639
9. Dokumentation, Varian, Inc.

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Bibliographie

1. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
2. Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc; 1984.
3. Hoffman B, Lefkowitz R. Catecholamines and Sympathomimetic Drugs. In: Gilman A, Rall T, Neis A, Taylor P, eds. *Goodman and Gillman's The Pharmacological Basis of Therapeutics*. 8th ed. New York, NY: McGraw Hill Pub Co; 1993:187-220.
4. Rowland M. Amphetamine blood and urine levels in man. *J Pharm Sci*. 1969;58(4):508-509..
5. Dimijian GG. Contemporary drug abuse. In: Goth A. ed. *Medical Pharmacology: Principles and Concepts*. 11th ed. St. Louis, Mo: CV Mosby Company; 1984:336-365..
6. Dring LG, Smith RL, Williams RT. The metabolic fate of amphetamine in man and other species. *Biochem J*. 1970;116:425-435.
7. Caldwell J. The metabolism of amphetamines in mammals. *Drug Metab Rev*. 1976;5(2):219-280.
8. Beckett AH, Rowland M. Urinary excretion kinetics of amphetamine in man. *J Pharm Pharmacol*. 1965;17:628-639
9. Documentation de Varian, Inc.

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- Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
- Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc: 1984.
- Hoffman B, Lefkowitz R. Catecholamines and Sympathomimetic Drugs. In: Gilman A, Rall T, Neis A, Taylor P, eds. *Goodman and Gillman's The Pharmacological Basis of Therapeutics*. 8th ed. New York, NY: McGraw Hill Pub Co; 1993:187-220.
- Rowland M. Amphetamine blood and urine levels in man. *J Pharm Sci*. 1969;58(4):508-509.
- Dimijian GG. Contemporary drug abuse. In: Goth A. ed. *Medical Pharmacology: Principles and Concepts*. 11th ed. St. Louis, Mo: CV Mosby Company; 1984:336-365..
- Dring LG, Smith RL, Williams RT. The metabolic fate of amphetamine in man and other species. *Biochem J*. 1970;116:425-435.
- Caldwell J. The metabolism of amphetamines in mammals. *Drug Metab Rev*. 1976;5(2):219-280.
- Beckett AH, Rowland M. Urinary excretion kinetics of amphetamine in man. *J Pharm Pharmacol*. 1965;17:628-639
- Datos de archivo, Varian, Inc.

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Riferimenti

- Hawks RL, Chiang CN (a cura di). Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
- Blum K. *Handbook of Abusable Drugs*. Prima ed. New York, NY: Gardner Press Inc: 1984.
- Hoffman B, Lefkowitz R. Catecholamines and Sympathomimetic Drugs. In: Gilman A, Rall T, Neis A, Taylor P (a cura di) *Goodman and Gillman's The Pharmacological Basis of Therapeutics*. ottava ed. New York, NY: McGraw Hill Pub Co; 1993:187-220.
- Rowland M. Amphetamine blood and urine levels in man. *J Pharm Sci*. 1969;58(4):508-509..
- Dimijian GG. Contemporary drug abuse. In: Goth A. (a cura di) *Medical Pharmacology: Principles and Concepts*. undicesima ed. St. Louis, Mo: CV Mosby Company; 1984:336-365..
- Dring LG, Smith RL, Williams RT. The metabolic fate of amphetamine in man and other species. *Biochem J*. 1970;116:425-435.
- Caldwell J. The metabolism of amphetamines in mammals. *Drug Metab Rev*. 1976;5(2):219-280.
- Beckett AH, Rowland M. Urinary excretion kinetics of amphetamine in man. *J Pharm Pharmacol*. 1965;17:628-639
- Dati presenti nel file, Varian, Inc.

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Bibliografia

1. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
2. Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc: 1984.
3. Hoffman B, Lefkowitz R. Catecholamines and Sympathomimetic Drugs. In: Gilman A, Rall T, Neis A, Taylor P, eds. *Goodman and Gilman's The Pharmacological Basis of Therapeutics*. 8th ed. New York, NY: McGraw Hill Pub Co; 1993:187-220.
4. Rowland M. Amphetamine blood and urine levels in man. *J Pharm Sci*. 1969;58(4):508-509..
5. Dimijian GG. Contemporary drug abuse. In: Goth A. ed. *Medical Pharmacology: Principles and Concepts*. 11th ed. St. Louis, Mo: CV Mosby Company; 1984:336-365..
6. Dring LG, Smith RL, Williams RT. The metabolic fate of amphetamine in man and other species. *Biochem J*. 1970;116:425-435.
7. Caldwell J. The metabolism of amphetamines in mammals. *Drug Metab Rev*. 1976;5(2):219-280.
8. Beckett AH, Rowland M. Urinary excretion kinetics of amphetamine in man. *J Pharm Pharmacol*. 1965;17:628-639
9. Ficheiros da Varian, Inc.

www.varian-onsite.com

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Dentro dos E.U.A.: 1-800-737-9667

Fora dos E.U.A.: 00 1 (949) 770-9381, extensão 2

DAT.TechSupport@varianinc.com

 Varian, Inc.
Lake Forest, CA 92630 USA

 Varian B.V.
Herculesweg 8, 4338 PL Middelburg
 Holanda

Número da patente: 5,770,458

16148 1

© 2003, Varian, Inc.

VARIAN