English

Intended Use

OnTrak TesTstik™ for Barbiturates is an *in vitro* diagnostic test intended for professional use for the qualitative detection of barbiturates in urine at or above a cutoff concentration of 200 ng/mL. OnTrak TesTstik for Barbiturates is not intended for over-the-counter sale.

OnTrak TesTstik provides only a preliminary analytical test result.

A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmation method.¹

Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Summary

The barbiturates, a class of drugs derived from barbituric acid, (malonylurea), are sedative hypnotics with central nervous system (CNS)-depressant activity. As CNS depressants, they are classified relative to their durations of action (ultra short, short, intermediate, and long-acting).² Barbiturates are used medically as sedatives to reduce emotional tension and induce sleep, and in certain types of epilepsy, they reduce seizure frequency by raising the seizure threshold.³ Excessive dosages may cause impaired motor coordination (slurred speech, loss of balance), perceptual alterations, faulty judgement, inflated perceptions of performance⁴ and disinhibition euphoria.^{3,5} Overdoses can result in stupor, coma, and death. The combined use of the barbiturates with alcohol, opiates, or other CNS depressants can result in fatal, additive respiratory depression.³ Oral administration is most common, although the barbiturates may be injected intravenously or intramuscularly. Following ingestion, they are rapidly absorbed from the stomach and enter the circulation. Their resulting distribution and concentration in various tissues is largely dependent on the lipid solubility and protein-binding characteristics of the different barbiturates.^{3,5,6} Most of the barbiturates are metabolized by the liver via side-chain oxidation, nitrogen dealkylation, and/or desulfuration of thiobarbiturates. The extent of liver metabolism is drug-dependent; secobarbital is extensively oxidized to a series of pharmacologically inactive metabolites, while phenobarbital and barbital are essentially excreted unchanged in the urine. The barbiturates are excreted as active drug/metabolite mixes whose ratios and concentrations depend on the specific barbiturate in question.⁷

Test Principle

OnTrak TesTstik assays are based on the principles of microparticle capture inhibition. The test relies on the competition between drug, which may be present in the urine being tested, and drug conjugate immobilized on membrane for binding to antibody coated onto colored microparticles. When the TesTstik is immersed in the urine specimen, some of the specimen is absorbed into the TesTstik sample pad. The absorbed specimen travels through a reagent strip contained in the device by capillary action. In the reagent strip, the specimen rehydrates and mobilizes antibody-coated blue microparticles. The microparticle-urine suspension continues to migrate through the reagent strip and comes in contact with the immobilized drug conjugate. In the absence of drug in the urine, the antibody-coated microparticles bind to the drug conjugate and a blue band is formed at the result window (“negative” sign).

When drug is present in the specimen, it binds to the antibody-coated microparticles. If sufficient drug is present, the microparticles are inhibited from binding the drug conjugate and no blue band is formed at the result window. A presumptive positive (“non-negative”) sample causes the membrane to remain white (“positive” sign).

An additional antibody/antigen reaction occurs at the “TEST VALID” area. The “TEST VALID” blue band forms when anti-bovine serum albumin (anti-BSA) antibodies, which are imbedded in the reagent membrane, bind to the BSA antigen on the blue microparticles. The presence of the “TEST VALID” band indicates that the test has completed, the reagents are viable, and the results are ready to interpret.

Reagents

1. Blue-dyed microparticles coated with mouse monoclonal anti-barbiturate analogue in a buffered solution containing preservative and dried onto a membrane.
2. Drug conjugates immobilized on membrane.
3. Mouse monoclonal anti-BSA antibody immobilized on membrane.

Precautions and Warnings

For *in vitro* diagnostic use. Exercise the normal precautions for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines.

Storage and Stability

OnTrak TesTstik should be stored at room temperature, 15 - 30 °C (59 - 86 °F), in the original foil pouch. If the foil pouch containing the TesTstik is damaged (e.g., a hole or tear), do not use that TesTstik. Do not use the TesTstik after the expiration date printed on the foil pouch.

Specimen Collection and Preparation

OnTrak TesTstik is formulated for use with urine specimens. Fresh urine specimens do not require any special handling or pretreatment. No additives or preservatives are required. It is recommended that specimens be stored refrigerated (2 - 8 °C or 36 - 46 °F) and tested within two days of collection. For prolonged storage, freezing is recommended.

Specimens containing human-sourced materials should be handled as if potentially infectious. Use safe laboratory procedures such as those outlined in *Biosafety in Microbiological and Biomedical Laboratories* (HHS Publication Number [CDC] 93-8395).

Testing Procedure

Materials provided

50 individually packaged OnTrak TesTstik for Barbiturates sticks

Ancillary materials

	Part No.	Quantity
Evidence Tape	42212	100 strips/pkg.
Latex Gloves	42210	50 pairs
Collection Cups	47123	100 cups/box
Adulteration Strips	3111407	25 strips/box

Contact Customer Technical Support for information regarding commercially available control material.

Assay

OnTrak TesTstik is ready for use directly from the foil pouch and must be used within eight hours after removal from the pouch. This time frame may be shorter in high humidity climates. Do not open the foil pouch unless it is at room temperature (15 - 30 °C or 59 - 86 °F).

Specimens should be between room temperature and body temperature 15 - 37 °C (59 - 99 °F). If the specimen was refrigerated, allow it to reach room temperature before testing. To ensure a sufficient quantity of sample for confirmation testing (as necessary), collect a minimum of 30 mL of sample.

- To expose the sample pad, pull the protective slide back until it stops (up to the "Dip" line).
- Gently immerse the exposed sample pad into the specimen until the specimen reaches the dip line. Do not immerse the entire device into the specimen. Keep the OnTrak TesTstik for Barbiturates in the specimen for approximately 10 seconds.
- Remove the OnTrak TesTstik for Barbiturates from the specimen and while holding the OnTrak TesTstik for Barbiturates over the specimen, push the protective slide back over the sample pad until it stops. Allow the test to proceed until a blue band is observed in the "TEST VALID" window (see Results section below). The OnTrak TesTstik for Barbiturates may be placed down or held while the test is developing. There is no required timing step. Note: If the Test Valid band does not appear, re-dip the OnTrak TesTstik for Barbiturates for approximately 10 seconds. Allow to sit undisturbed for an additional 1-3 minutes. Once the Test Valid band forms, proceed to step 4.
- Remove the tab that covers the result window from the protective slide. Read and record the result.

Quality Control

Quality control testing at regular intervals is good laboratory practice and may be required by federal, state, or local guidelines. Always check with the appropriate licensing or accrediting bodies to ensure your quality control program meets the established standards.

Internal control: Each OnTrak TesTstik BAR has built-in process controls. The "TEST VALID" band is an internal POSITIVE process control. A blue "TEST VALID" band should always appear if an adequate sample volume is used, the sample and reagent are wicking on the membrane, and the reagents in the "TEST VALID" band and the conjugate-color indicator are working. In addition, the background in the result window should become clear and provide a distinct result. If smears of color appear in the background, nonspecific binding of the antibody may be occurring, making the test result difficult to interpret. The clearing of the result window background may be considered an internal NEGATIVE process control. If the "TEST VALID" band does not appear in the test window, or the background does not clear, then the test is not valid and a new test must be performed.

External control: The use of external positive and negative controls is recommended to test each shipment of product or when a new lot is used. Contact Customer Technical Support for information regarding commercially available control material. The controls should produce the expected results. If expected results are not obtained, do not run test specimens, and contact Customer Technical Support. Follow the appropriate federal, state and local guidelines when running external controls.

Results

Once a full, partial or broken blue band is formed in the "TEST VALID" window (in approximately 3 minutes) remove the cover tab and interpret results as either presumptive positive or negative. Results may be interpreted for up to 30 minutes after the "TEST VALID" band forms.

A negative result is the presence of a blue band ("negative" sign) anywhere in the test result window. Since a negative sample may give a faint or incomplete color band, **any** full, partial or broken band in the test result window indicates a negative result.

A presumptive positive result (drug present at or above the cutoff) is the lack of a blue band in the test result window. The window appears as a white plus (+) sign.

Limitations

See Specific Performance Characteristics for information on substances tested for cross-reactivity in this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors). Adulteration of the sample with oxidizing agents and/or acidic or alkaline agents can cause erroneous results. If adulteration is suspected, another sample should be collected and a new OnTrak TesTstik used. For information on how adulterants may affect the performance of this assay, contact Customer Technical Support.

A presumptive positive result with this assay indicates the presence of drug or drug metabolites in urine but does not reflect the degree of intoxication.

Specific Performance Characteristics⁸

Accuracy

OnTrak TesTstik for Barbiturates was evaluated using specimens screened by an automated immunoassay and confirmed positive by GC/MS (at a 200 ng/mL cutoff). Fifty samples positive for barbiturates were positive by OnTrak TesTstik (100%).

		GC/MS	
		+	
TesTstik	+	50	
	-	0	

One hundred urine samples, obtained from a clinical laboratory and screened negative by an automated immunoassay relative to a 200 ng/mL cutoff for barbiturates were evaluated and found negative using OnTrak TesTstik (100%).

All positive and negative samples were also assayed by, and compared to, Abuscreen OnTrak for Barbiturates (200 ng/mL). All samples demonstrated 100% agreement between the two assays.

		Abuscreen OnTrak	
		+	-
TesTstik	+	50	0
	-	0	100

Precision

The precision of OnTrak TesTstik for Barbiturates was determined by testing various concentrations of drug standards for the primary drugs or drug metabolites detected by this system. All lots are required to minimally perform at a greater than 95% confidence level that a negative result will be attained with drug at 25% (¼ x) of its cutoff concentration, and a greater than 95% confidence level that a positive result will be attained with drug at 150% (1½ x) of its cutoff concentration.

Specificity

The following structurally similar compounds for barbiturates were tested for cross-reactivity with OnTrak TesTstik for Barbiturates. Each compound tested was prepared in normal human urine. The results are expressed as that amount of the compound capable of giving a result equivalent to the cutoff.

Barbiturates-related compounds	Approx. ng/mL equivalent to 200 ng/mL secobarbital	Approx. percent cross-reactivity
Allobarbital	200	100
Aprobarbital	200	100
Butalbital	200	100
Cyclopentobarbital	200	100
Butabarbital	400	50
Pentobarbital	400	50
Amobarbital	750	27
Barbital	1000	20
p-Hydroxyphenobarbital	1000	20
Phenobarbital	1000	20
Barbituric acid	>100,000	<0.2
1,3-Dimethylbarbituric acid	>100,000	<0.2
Hexobarbital	>100,000	<0.2
Mephobarbital	>100,000	<0.2

Cross-Reactivity with Additional Drugs

The following compounds were tested at 100,000 ng/mL and found not to cross-react with the OnTrak TesTstik for Barbiturates assay.

Acetaminophen	Benzoylcegonine
Acetylsalicylic acid	Brompheniramine
Aminopyrine	Caffeine
Amisriptyline	Calcium hypochlorite
Amobarbital	Chlordiazepoxide
Amphetamine	Chloroquine
Ampicillin	Chlorpheniramine
Ascorbic acid	Chlorpromazine
Aspartame	Clemastine
Atropine	Cocaine
Benzocaine	Codeine

1 Gently tear open the foil pouch to remove the stick. Pull the protective slide back to expose the sample pad.

2 Immerse the sample pad into the specimen until the specimen reaches the DIP LINE. Hold the stick in the specimen for approximately 10 seconds.

3 Remove the stick from the specimen. While holding the stick over the specimen cup, push the protective slide back over the sample pad.


Test complete—
ready to interpret
results

Test NOT complete—
DO NOT interpret
results

4 Allow the test to proceed until a blue band forms in the “TEST VALID” window (see example on the LEFT). Test is NOT complete until a blue band forms in the “TEST VALID” window (see example on the RIGHT). Once the “TEST VALID” line appears, remove the cover tab and interpret results.

5 A blue band in the result window indicates a NEGATIVE result.

6 No color in the result window indicates a POSITIVE result. The window appears as a white plus (+) sign.

Cyclizine
Desipramine
Dextromethorphan
Dextropropoxyphene
Diazepam
Diphenhydramine
Diphenylhydantoin
Dopamine
Doxylamine
Ecgonine
Ecgonine methyl ester
Ephedrine
Epinephrine
Erythromycin
Estriol
Ethinylestradiol
Fenopropfen
Furosemide
Gentisic acid
Glutethimide
Guaiacol glycerol ether
Hydrochlorothiazide
p-Hydroxyamphetamine
3-Hydroxytyramine
Ibuprofen
Imipramine
Isoproterenol
Ketamine
LSD
Lidocaine
MDA
MDMA
Melanin
Meperidine
Methadone
Methamphetamine
Methapyrilene
Methaqualone
Methylphenidate

Methyprylon
Morphine Sulfate
Naloxone
Naltrexone
Naproxen
Niacinamide
Nordiazepam
Norethindrone
Norpseudoephedrine
Nortriptyline
Oxazepam
Penicillin G
Phencyclidine
 β -Phenethylamine
Phenothiazine
Phentermine
Phenylbutazone
Phenylephrine
Phenylpropanolamine
Phenyltoloxamine
Procaine
Procyclidine
Promethazine
Pseudoephedrine
Quinidine
Quinine
Sulindac
 Δ^9 -THC-9-carboxylic acid
Tetracycline
Tetrahydrozoline
Thiordazine
Trifluoperazine
Trihexyphenidyl
Trimipramine
Tripelennamine
Tyramine
Verapamil
Zomepirac

References

1. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73. 1986.
2. Clark WG, Brater DC, Johnson AR, eds. *Goth's Medical Pharmacology*. 12th ed. St. Louis: CV Mosby Company; 1988.
3. Wesson DR, Smith DE. *Barbiturates: Their Use, Misuse, and Abuse*. New York: Human Sciences Press; 1977.
4. Smith GM, Beecher HK. Amphetamine, secobarbital, and athletic performance: III. Quantitative effects on judgement. *JAMA*. 1960;172:1623-1629.
5. Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc; 1984.
6. Robinson AE, McDowall RD. The distribution of amylobarbitone, butobarbitone, pentobarbitone, and quinalbarbitone, and the hydroxylated metabolites in man. *J Pharm Pharmacol*. 1979;31:357-365.
7. Mark LC. Metabolism of barbiturates in man. *Clin Pharmacol Ther*. 1963;4: 504-530.
8. Data on file at Varian, Inc.

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Literatur

- Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73. 1986.
- Clark WG, Brater DC, Johnson AR, eds. Goth's Medical Pharmacology. 12th ed. St. Louis: CV Mosby Company; 1988.
- Wesson DR, Smith DE. Barbiturates: Their Use, Misuse, and Abuse. New York: Human Sciences Press; 1977.
- Smith GM, Beecher HK. Amphetamine, secobarbital, and athletic performance: III. Quantitative effects on judgement. *JAMA*. 1960;172:1623-1629.
- Blum K. Handbook of Abusable Drugs. 1st ed. New York, NY: Gardner Press Inc; 1984.
- Robinson AE, McDowall RD. The distribution of amylobarbitone, butobarbitone, pentobarbitone, and quinalbarbitone, and the hydroxylated metabolites in man. *J Pharm Pharmacol*. 1979;31:357-365.
- Mark LC. Metabolism of barbiturates in man. *Clin Pharmacol Ther*. 1963;4: 504-530.
- Dokumentation, Varian, Inc.

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Bibliographie

1. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73. 1986.
2. Clark WG, Brater DC, Johnson AR, eds. *Goth's Medical Pharmacology*. 12th ed. St. Louis: CV Mosby Company; 1988.
3. Wesson DR, Smith DE. Barbiturates: Their Use, Misuse, and Abuse. New York: Human Sciences Press; 1977.
4. Smith GM, Beecher HK. Amphetamine, secobarbital, and athletic performance: III. Quantitative effects on judgement. *JAMA*. 1960;172:1623-1629.
5. Blum K. Handbook of Abusable Drugs. 1st ed. New York, NY: Gardner Press Inc; 1984.
6. Robinson AE, McDowall RD. The distribution of amylobarbitone, butobarbitone, pentobarbitone, and quinalbarbitone, and the hydroxylated metabolites in man. *J Pharm Pharmacol*. 1979;31:357-365.
7. Mark LC. Metabolism of barbiturates in man. *Clin Pharmacol Ther*. 1963;4: 504-530.
8. Documentation de Varian, Inc.

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- Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73. 1986.
- Clark WG, Brater DC, Johnson AR, eds. Goth's Medical Pharmacology. 12th ed. St. Louis: CV Mosby Company; 1988.
- Wesson DR, Smith DE. Barbiturates: Their Use, Misuse, and Abuse. New York: Human Sciences Press; 1977.
- Smith GM, Beecher HK. Amphetamine, secobarbital, and athletic performance: III. Quantitative effects on judgement. *JAMA*. 1960;172:1623-1629.
- Blum K. Handbook of Abusable Drugs. 1st ed. New York, NY: Gardner Press Inc; 1984.
- Robinson AE, McDowall RD. The distribution of amylobarbitone, butobarbitone, pentobarbitone, and quinalbarbitone, and the hydroxylated metabolites in man. *J Pharm Pharmacol*. 1979; 31:357-365.
- Mark LC. Metabolism of barbiturates in man. *Clin Pharmacol Ther*. 1963;4: 504-530.
- Datos de archivo, Varian, Inc.

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Riferimenti

1. Hawks RL, Chiang CN (a cura di) Urine testing for drugs of abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73. 1986.
2. Clark WG, Brater DC, Johnson AR (a cura di). Goth's Medical Pharmacology. dodicesima ed. St. Louis, MO: CV Mosby Company; 1988.
3. Wesson DR, Smith DE. Barbiturates: Their Use, Misuse, and Abuse. New York: Human Sciences Press; 1977.
4. Smith GM, Beecher HK. Amphetamine, secobarbital, and athletic performance: III. Quantitative effects on judgement. *JAMA*. 1960;172:1623-1629.
5. Blum K. Handbook of Abusable Drugs. prima ed. New York, NY: Gardner Press Inc; 1984.
6. Robinson AE, McDowall RD. The distribution of amylobarbitone, butobarbitone, pentobarbitone, and quinalbarbitone, and the hydroxylated metabolites in man. *J Pharm Pharmacol*. 1979;31: 357-365.
7. Mark LC. Metabolism of barbiturates in man. *Clin Pharmacol Ther*. 1963;4: 504-530.
8. Dati presenti nel file, Varian, Inc.

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Bibliografia

1. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73. 1986.
2. Clark WG, Brater DC, Johnson AR, eds. Goth's Medical Pharmacology. 12th ed. St. Louis, MO: CV Mosby Company; 1988.
3. Wesson DR, Smith DE. Barbiturates: Their Use, Misuse, and Abuse. New York: Human Sciences Press; 1977.
4. Smith GM, Beecher HK. Amphetamine, secobarbital, and athletic performance: III. Quantitative effects on judgement. *JAMA*. 1960;172:1623-1629.
5. Blum K. Handbook of Abusable Drugs. 1st ed. New York, NY: Gardner Press Inc; 1984.
6. Robinson AE, McDowall RD. The distribution of amylobarbitone, butobarbitone, pentobarbitone, and quinalbarbitone, and the hydroxylated metabolites in man. *J Pharm Pharmacol*. 1979;31:357-365.
7. Mark LC. Metabolism of barbiturates in man. *Clin Pharmacol Ther*. 1963;4: 504-530.
8. Ficheiros da Varian, Inc.

www.varian-onsite.com

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Dentro dos E.U.A.: 1-800-737-9667

Fora dos E.U.A.: 00 1 (949) 770-9381, extensão 2

DAT.TechSupport@varianinc.com

 Varian, Inc.
Lake Forest, CA 92630 USA

 Varian B.V.
Herculesweg 8, 4338 PL Middelburg
 Holanda

Número da patente: 5,770,458

16149 1

© 2003, Varian, Inc.

VARIAN