English

Intended Use

OnTrak TesTstik™ for Cocaine is an *in vitro* diagnostic test intended for professional use for the qualitative detection of cocaine metabolite in urine at or above a cutoff concentration of 300 ng/mL. OnTrak TesTstik for Cocaine is not intended for over-the-counter sale.

OnTrak TesTstik provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmation method.¹ Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Summary

Cocaine, a natural product found in the leaves of the coca plant, is a potent central nervous system (CNS) stimulant and a local anesthetic. Its pharmacological effects are identical to those of the amphetamines (also CNS stimulants), though cocaine has a shorter duration of action.² Cocaine induces euphoria, confidence and a sense of increased energy in the user; these psychological effects are accompanied by increased heart rate, dilation of pupils, fever, tremors, and sweating. The “crash” following a cocaine “high” is profound, ranging from irritability, lassitude, and the desire for more drug, to anxiety, hallucinations, and paranoia.^{3,4} Users may resort to other drugs at this time to relieve the depressive effects of the “crash.”²

Cocaine is traditionally administered intranasally or smoked in its purer, free-base form. Oral ingestion is ineffective, as cocaine is broken down in the gastrointestinal tract. It is absorbed readily across the mucous membranes of the nose and lungs into the circulation. Its effects are intense but short-lived. Cocaine is rapidly inactivated by hydrolysis of its ester linkages.^{1,5,6} Blood cholinesterases hydrolyze cocaine to ecgonine methyl ester, while hydrolysis of the parent drug to benzoylecgonine is thought to be nonenzymatic. Both of these metabolites may be further hydrolyzed to ecgonine. Unmetabolized cocaine has an affinity for fatty tissue and rapidly enters the brain. Cocaine metabolites, however, are more water soluble and are readily excreted in the urine along with some portion of unchanged drug.^{5,7} The prominent benzoylecgonine metabolite is the primary urinary marker for detecting cocaine use.^{1,5}

Tolerance has been observed with some chronic, high-dose users. Physical dependence does not appear to occur in abusers, although the development of strong psychological dependence is well known. Cessation of drug use may result in depression, hallucinations, and in extreme cases, psychosis.²

Test Principle

OnTrak TesTstik assays are based on the principles of microparticle capture inhibition. The test relies on the competition between drug, which may be present in the urine being tested, and drug conjugate immobilized on membrane for binding to antibody coated onto colored microparticles. When the TesTstik is immersed in the urine specimen, some of the specimen is absorbed into the TesTstik sample pad. The absorbed specimen travels through a reagent strip contained in the device by capillary action. In the reagent strip, the specimen rehydrates and mobilizes antibody-coated blue microparticles. The microparticle-urine suspension continues to migrate through the reagent strip and comes in contact with the immobilized drug conjugate. In the absence of drug in the urine, the antibody-coated microparticles bind to the drug conjugate and a blue band is formed at the result window (“negative” sign).

When drug is present in the specimen, it binds to the antibody-coated microparticles. If sufficient drug is present, the microparticles are inhibited from binding the drug conjugate and no blue band is formed at the result window. A presumptive positive (“non-negative”) sample causes the membrane to remain white (“positive” sign).

An additional antibody/antigen reaction occurs at the “TEST VALID” area. The “TEST VALID” blue band forms when anti-bovine serum albumin (anti-BSA) antibodies, which are imbedded in the reagent membrane, bind to the BSA antigen on the blue microparticles. The presence of the “TEST VALID” band indicates that the test has completed, the reagents are viable, and the results are ready to interpret.

Reagents

1. Blue-dyed microparticles coated with mouse monoclonal anti-benzoylecgonine antibody and BSA in a buffered solution containing preservative and dried onto a membrane.
2. Drug conjugates immobilized on membrane.
3. Mouse monoclonal anti-BSA antibody immobilized on membrane.

Precautions and Warnings

For *in vitro* diagnostic use. Exercise the normal precautions for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines.

Storage and Stability

OnTrak TesTstik should be stored at room temperature, 15 - 30 °C (59 - 86 °F), in the original foil pouch. If the foil pouch containing the TesTstik is damaged (e.g., a hole or tear), do not use that TesTstik. Do not use the TesTstik after the expiration date printed on the foil pouch.

Specimen Collection and Preparation

OnTrak TesTstik is formulated for use with urine specimens. Fresh urine specimens do not require any special handling or pretreatment. No additives or preservatives are required. It is recommended that specimens be stored refrigerated (2 - 8 °C or 36 - 46 °F) and tested within two days of collection. For prolonged storage, freezing is recommended.

Specimens containing human-sourced materials should be handled as if potentially infectious. Use safe laboratory procedures such as those outlined in *Biosafety in Microbiological and Biomedical Laboratories* (HHS Publication Number [CDC] 93-8395).

Testing Procedure

Materials provided

50 individually packaged OnTrak TesTstik for Cocaine sticks

Ancillary materials

	Part No.	Quantity
Evidence Tape	42212	100 strips/pkg.
Latex Gloves	42210	50 pairs
Collection Cups	47123	100 cups/box
Adulteration Strips	3111407	25 strips/box

Contact Customer Technical Support for information regarding commercially available control material.

Assay

OnTrak TesTstik is ready for use directly from the foil pouch and must be used within eight hours after removal from the pouch. This time frame may be shorter in high humidity climates. Do not open the foil pouch unless it is at room temperature (15 - 30 °C or 59 - 86 °F).

Specimens should be between room temperature and body temperature 15 - 37 °C (59 - 99 °F). If the specimen was refrigerated, allow it to reach room temperature before testing. To ensure a sufficient quantity of sample for confirmation testing (as necessary), collect a minimum of 30 mL of sample.

1. To expose the sample pad, pull the protective slide back until it stops (up to the "Dip" line).
2. Gently immerse the exposed sample pad into the specimen until the specimen reaches the dip line. Do not immerse the entire device into the specimen. Keep the OnTrak TesTstik for Cocaine in the specimen for approximately 10 seconds.
3. Remove the OnTrak TesTstik for Cocaine from the specimen and while holding the OnTrak TesTstik for Cocaine over the specimen, push the protective slide back over the sample pad until it stops. Allow the test to proceed until a blue band is observed in the "TEST VALID" window (see Results section below). The OnTrak TesTstik for Cocaine may be placed down or held while the test is developing. There is no required timing step. Note: If the Test Valid band does not appear, re-dip the OnTrak TesTstik for Cocaine for approximately 10 seconds. Allow to sit undisturbed for an additional 1-3 minutes. Once the Test Valid band forms, proceed to step 4.
4. Remove the tab that covers the result window from the protective slide. Read and record the result.

Quality Control

Quality control testing at regular intervals is good laboratory practice and may be required by federal, state, or local guidelines. Always check with the appropriate licensing or accrediting bodies to ensure your quality control program meets the established standards.

Internal control: Each OnTrak TesTstik COC has built-in process controls. The "TEST VALID" band is an internal POSITIVE process control. A blue "TEST VALID" band should always appear if an adequate sample volume is used, the sample and reagent are wicking on the membrane, and the reagents in the "TEST VALID" band and the conjugate-color indicator are working. In addition, the background in the result window should become clear and provide a distinct result. If smears of color appear in the background, nonspecific binding of the antibody may be occurring, making the test result difficult to interpret. The clearing of the result window background may be considered an internal NEGATIVE process control. If the "TEST VALID" band does not appear in the test window, or the background does not clear, then the test is not valid and a new test must be performed.

External control: The use of external positive and negative controls is recommended to test each shipment of product or when a new lot is used. Contact Customer Technical Support for information regarding commercially available control material. The controls should produce the expected results. If expected results are not obtained, do not run test specimens, and contact Customer Technical Support. Follow the appropriate federal, state and local guidelines when running external controls.

Results

Once a full, partial or broken blue band is formed in the "TEST VALID" window (in approximately 3 minutes) remove the cover tab and interpret results as either presumptive positive or negative. Results may be interpreted for up to 30 minutes after the "TEST VALID" band forms.

A negative result is the presence of a blue band ("negative" sign) anywhere in the test result window. Since a negative sample may give a faint or incomplete color band, **any** full, partial or broken band in the test result window indicates a negative result.

A presumptive positive result (drug present at or above the cutoff) is the lack of a blue band in the test result window. The window appears as a white plus (+) sign.

Limitations

See Specific Performance Characteristics for information on substances tested for cross-reactivity in this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors). Adulteration of the sample with oxidizing agents and/or acidic or alkaline agents can cause erroneous results. If adulteration is suspected, another sample should be collected and a new OnTrak TesTstik used. For information on how adulterants may affect the performance of this assay, contact Customer Technical Support.

A presumptive positive result with this assay indicates the presence of drug or drug metabolites in urine but does not reflect the degree of intoxication.

Specific Performance Characteristics⁸

Accuracy

OnTrak TesTstik for Cocaine was evaluated using specimens screened by an automated immunoassay and confirmed positive by GC/MS (at a 150 ng/mL cutoff). Fifty samples positive for cocaine were positive by OnTrak TesTstik (100%).

		GC/MS
		+
TesTstik	+	50
	-	0

One hundred six urine samples, obtained from a clinical laboratory and screened negative by an automated immunoassay relative to a 300 ng/mL cutoff for cocaine were evaluated and found negative using OnTrak TesTstik.

All positive and negative samples were also assayed by, and compared to, Abuscreen OnTrak for Cocaine (300 ng/mL). All samples demonstrated 100% agreement between the two assays.

		Abuscreen OnTrak	
		+	-
TesTstik	+	50	0
	-	0	106

Precision

The precision of OnTrak TesTstik for Cocaine was determined by testing various concentrations of drug standards for the primary drugs or drug metabolites detected by this system. All lots are required to minimally perform with a greater than 95% confidence level that negative results will be attained with drugs at 25% (¼ x) of their respective cutoff concentrations, and a greater than 95% confidence level that positive results will be attained with drugs at 150% (1½ x) of their respective cutoff concentrations.

Specificity

The following structurally similar compounds for cocaine were tested for cross-reactivity with OnTrak TesTstik for Cocaine. Each compound tested was prepared in normal human urine. The results are expressed as that amount of the compound capable of giving a result equivalent to the cutoff of the assay.

Cocaine-related compounds	Approx. ng/mL equivalent to 300 ng/mL benzoyllecgonine	Approx. percent cross-reactivity
Cocaine HCl	40,000	0.75
Ecgonine HCl	>100,000	<0.3
Ecgonine methyl ester HCl	>100,000	<0.3

1 Gently tear open the foil pouch to remove the stick. Pull the protective slide back to expose the sample pad.

2 Immerse the sample pad into the specimen until the specimen reaches the DIP LINE. Hold the stick in the specimen for approximately 10 seconds.

3 Remove the stick from the specimen. While holding the stick over the specimen cup, push the protective slide back over the sample pad.


Test complete— ready to interpret results

Test NOT complete— DO NOT interpret results

4 Allow the test to proceed until a blue band forms in the “TEST VALID” window (see example on the LEFT). Test is NOT complete until a blue band forms in the “TEST VALID” window (see example on the RIGHT). Once the “TEST VALID” line appears, remove the cover tab and interpret results.

5 A blue band in the result window indicates a NEGATIVE result.

6 No color in the result window indicates a POSITIVE result. The window appears as a white plus (+) sign.

Cross-Reactivity with Additional Drugs

The following compounds were tested at 100,000 ng/mL and found not to cross-react with the OnTrak TesTstik for Cocaine assay.

Acetaminophen	Estriol
Acetylsalicylic acid	Fenoprofen
Aminopyrine	Furosemide
Amitriptyline	Gentisic acid
Amobarbital	Glutethimide
Amphetamine	Guaiacol glycerol ether
Ampicillin	Hydrochlorothiazide
Ascorbic acid	<i>p</i> -Hydroxyamphetamine
Aspartame	Ibuprofen
Atropine	Imipramine
Benzocaine	Isoproterenol
Benzphetamine	Ketamine
Butabarbital	LSD
Caffeine	Lidocaine
Calcium hypochlorite	MDA
Chlordiazepoxide	MDMA
Chloroquine	Melanin
Chlorpheniramine	Meperidine
Chlorpromazine	Methadone
Codeine	Methamphetamine
Dextromethorphan	Methaqualone
Dextropropoxyphene	Methylphenidate
Diazepam	Methyprylon
Diphenhydramine	Morphine Sulfate
Diphenylhydantoin	Naloxone
Dopamine	Naltrexone
Doxylamine	Naproxen
Ephedrine	Niacinamide
Epinephrine	Nordiazepam
Erythromycin	Norethindrone

Norpseudoephedrine	Promethazine
Oxazepam	Pseudoephedrine
Penicillin G	Quinidine
Pentobarbital	Quinine
Phencyclidine	Secobarbital
β -Phenethylamine	Sulindac
Phenobarbital	Δ^9 -THC-9-carboxylic acid
Phenothiazine	Tetracycline
Phentermine	Tetrahydrozoline
Phenylbutazone	Trifluoperazine
Phenylephrine	Tyramine
Phenylpropanolamine	Verapamil
Procaine	

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Riferimenti

1. Hawks RL, Chiang CN (a cura di) Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
2. Blum K. *Handbook of Abusable Drugs*. prima ed. New York, NY: Gardner Press Inc; 1984.
3. Wilson MC, et al. Acute pharmacological activity of intravenous cocaine in the rhesus monkey. *Psychopharmacol Comm.* 1976;2:251.
4. Rappolt RT, et al. Propranolol in the treatment of cardiopressor effects of cocaine. *N Engl J Med.* 1976;295:448.
5. Fish F, Wilson WDC. Excretion of cocaine and its metabolites in man. *J Pharm Pharmacol.* 1969;21:135S.
6. Hamilton HE, et al. Cocaine and benzoylecgonine excretion in humans. *J Forensic Sci.* 1977;22:697.
7. Stewart DJ, et al. Cocaine metabolism: Cocaine and norcocaine hydrolysis by liver and serum esterases. *Clin Pharmacol Ther.* 1979;25:464.
8. Dati presenti nel file, Varian, Inc.

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© 2003, Varian, Inc.
2. Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc; 1984.
8. Ficheiros da Varian, Inc.
Varian B.V.