English **Intended Use**

OnTrak TesTstik™ for Methamphetamine is an *in vitro* diagnostic test intended for professional use for the qualitative detection of methamphetamine in urine at or above a cutoff concentration of 500 ng/mL. OnTrak TesTstik for Methamphetamine is not intended for over-the-counter sale.

OnTrak TesTstik provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmation method.¹ Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Summary

Methamphetamine, a sympathomimetic drug chemically related to amphetamine and structurally similar to designer drugs, e.g., "Ecstasy," has prominent effects on the central nervous system (CNS). Methamphetamine intake suppresses appetite and stimulates alertness, energy, euphoria, and causes a rise in blood pressure and cardiac arrhythmias.² In toxic doses, however, unpleasant CNS symptoms may occur such as anxiety, agitation, hallucinations, psychosis, seizures and possible death.^{2,3} Chronic methamphetamine abusers are extremely dangerous; they suffer from insomnia, paranoia, schizophrenia and exhibit frequent violent behaviors.³

Methamphetamine is a powder that is generally administered orally or intravenously. A more powerful crystallized form ("ICE") is smoked.⁴ Once ingested, it is readily absorbed in the stomach and rapidly enters the bloodstream. It is then metabolized by the liver by ring hydroxylation, N-demethylation, deamination and conjugation.⁴ Under normal conditions, up to 43% of a d-methamphetamine dosage is excreted unchanged in the urine in the first 24 hours and amphetamine, the primary metabolite, is present in small amounts (4 - 7%).⁵ In acidic urine, up to 76% is present as parent drug compared with 2% under alkaline conditions.⁵ Certain factors, such as individual metabolism, diet, dose, frequency of abuse, and time of ingestion may influence the presence of the parent drug and its metabolites.⁶

Test Principle

OnTrak TesTstik assays are based on the principles of microparticle capture inhibition. The test relies on the competition between drug, which may be present in the urine being tested, and drug conjugate immobilized on membrane for binding to antibody-coated blue microparticles.

When the TesTstik is immersed in the urine specimen, some of the specimen is absorbed into the TesTstik sample pad. The absorbed specimen travels through a reagent strip contained in the device by capillary action. In the reagent strip, the specimen rehydrates and mobilizes antibody-coated blue microparticles. The microparticle-urine suspension continues to migrate through the reagent strip and comes in contact with the immobilized drug conjugate. In the absence of drug in the urine, the antibody-coated microparticles bind to the drug conjugate and a blue band is formed at the result window ("negative" sign).

When drug is present in the specimen, it binds to the antibody-coated microparticles. If sufficient drug is present, the microparticles are inhibited from binding the drug conjugate and no blue band is formed at the result window. A presumptive positive ("non-negative") sample causes the membrane to remain white ("positive" sign).

An additional antibody/antigen reaction occurs at the "TEST VALID" area. The "TEST VALID" blue band forms when anti-bovine serum albumin (anti-BSA) antibodies imbedded in the reagent membrane bind

to the BSA antigen on the blue microparticles. The presence of the "TEST VALID" band indicates that the test has completed, the reagents are viable, and the results are ready to interpret.

Reagents

1. Blue-dyed microparticles coated with mouse monoclonal anti-methamphetamine antibody, and bovine serum albumin in a buffered solution containing preservative and dried onto a membrane.
2. Drug conjugates immobilized on membrane.
3. Mouse monoclonal anti-BSA antibody immobilized on membrane.

Precautions and Warnings

For *in vitro* diagnostic use. Exercise the normal precautions for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines.

Storage and Stability

OnTrak TesTstik should be stored at room temperature, 15 - 30 °C (59 - 86 °F), in the original foil pouch. If the foil pouch containing the TesTstik is damaged (e.g., a hole or tear), do not use that TesTstik. Do not use the TesTstik after the expiration date printed on the foil pouch.

Specimen Collection and Preparation

OnTrak TesTstik is formulated for use with urine specimens. Fresh urine specimens do not require any special handling or pretreatment. No additives or preservatives are required. It is recommended that specimens be stored refrigerated (2 - 8 °C or 36 - 46 °F) and tested within two days of collection. For prolonged storage, freezing is recommended.

Specimens containing human-sourced materials should be handled as if potentially infectious. Use safe laboratory procedures such as those outlined in *Biosafety in Microbiological and Biomedical Laboratories* (HHS Publication Number [CDC] 93-8395).

Testing Procedure

Materials provided

50 individually packaged OnTrak TesTstik for Methamphetamine sticks

Ancillary materials

	Part No.	Quantity
Evidence Tape	42212	100 strips/pkg.
Latex Gloves	42210	50 pairs
Collection Cups	47123	100 cups/box
Adulteration Strips	3111407	25 strips/box

Contact Customer Technical Support for information regarding commercially available control material.

Assay

OnTrak TesTstik is ready for use directly from the foil pouch and must be used within eight hours after removal from the pouch. This time frame may be shorter in high humidity climates. Do not open the foil pouch unless it is at room temperature (15 - 30 °C or 59 - 86 °F).

Specimens should be between room temperature and body temperature 15 - 37 °C (59 - 99 °F). If the specimen was refrigerated, allow it to reach room temperature before testing. To ensure a sufficient quantity of sample for confirmation testing (as necessary), collect a minimum of 30 mL of sample.

1. To expose the sample pad, pull the protective slide back until it stops (up to the "Dip" line).
2. Gently immerse the exposed sample pad into the specimen until the specimen reaches the dip line. Do not immerse the entire device into the specimen. Keep the OnTrak TesTstik for Methamphetamine in the specimen for approximately 10 seconds.

- Remove the OnTrak TesTstik for Methamphetamine from the specimen and while holding the OnTrak TesTstik for Methamphetamine over the specimen, push the protective slide back over the sample pad until it stops. Allow the test to proceed until a blue band is observed in the "TEST VALID" window (see Results section below). The OnTrak TesTstik for Methamphetamine may be placed down or held while the test is developing. There is no required timing step. Note: If the Test Valid band does not appear, re-dip the OnTrak TesTstik for Methamphetamine for approximately 10 seconds. Allow to sit undisturbed for an additional 1-3 minutes. Once the Test Valid band forms, proceed to step 4.
- Remove the tab that covers the result window from the protective slide. Read and record the result.

Quality Control

Quality control testing at regular intervals is good laboratory practice and may be required by federal, state, or local guidelines. Always check with the appropriate licensing or accrediting bodies to ensure your quality control program meets the established standards.

Internal control: Each OnTrak TesTstik for Methamphetamine has built-in process controls. The "TEST VALID" band is an internal POSITIVE process control. A blue "TEST VALID" band should always appear if an adequate sample volume is used, the sample and reagent are wicking on the membrane, and the reagents in the "TEST VALID" band and the conjugate-color indicator are working. In addition, the background in the result window should become clear and provide a distinct result. If smears of color appear in the background, nonspecific binding of the antibody may be occurring, making the test result difficult to interpret. The clearing of the result window background may be considered an internal NEGATIVE process control. If the "TEST VALID" band does not appear in the test window, or the background does not clear, then the test is not valid and a new test must be performed.

External control: The use of external positive and negative controls is recommended to test each shipment of product or when a new lot is used. Contact Customer Technical Support for information regarding commercially available control material. The controls should produce the expected results. If expected results are not obtained, do not run test specimens, and contact Customer Technical Support. Follow the appropriate federal, state and local guidelines when running external controls.

Results

Once a full, partial or broken blue band is formed in the "TEST VALID" window (in approximately 3 minutes) remove the cover tab and interpret results as either positive or negative. Results may be interpreted for up to 30 minutes after the "TEST VALID" band forms.

A negative result is the presence of a blue band ("negative" sign) anywhere in the test result window. Since negative specimens may give a faint or incomplete color band, **any** full, partial or broken band in the test result window indicates a negative result.

A presumptive positive result (drug present at or above the cutoff) is the lack of a blue band in the test result window. The window appears as a white plus (+) sign.

Limitations

See Specific Performance Characteristics for information on substances tested for cross-reactivity in this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors). Adulteration of the sample with oxidizing agents and/or acidic or alkaline agents can cause erroneous results. If adulteration is suspected, another sample should be collected and a new OnTrak TesTstik used. For information on how adulterants may affect the performance of this assay, contact Customer Technical Support.

A presumptive positive result with this assay indicates the presence of drug or drug metabolites in urine but does not reflect the degree of intoxication.

Specific Performance Characteristics⁷

Accuracy

OnTrak TesTstik for Methamphetamine was evaluated and compared to the Syva Emit II immunoassay using specimens prescreened by an automated immunoassay and confirmed positive by GC/MS. Fifty-five

specimens positive for methamphetamine were evaluated with OnTrak TesTstik for Methamphetamine. All specimens were positive by both OnTrak TesTstik for Methamphetamine and Syva Emit II.

One hundred twenty urine specimens were obtained from a clinical laboratory and screened negative by an automated immunoassay relative to a 500 ng/mL cutoff for methamphetamine. A portion (12%) of the negative specimens were also confirmed negative by GC/MS analysis for methamphetamine. Of the 120 negative specimens, all were negative by both OnTrak TesTstik for Methamphetamine and Syva Emit II.

OnTrak TesTstik was evaluated and compared to GC/MS using specimens selected to be within approximately 25% of the immunological cutoff (500 ng/mL). A total of 30 specimens were tested. Eight specimens were negative and 22 were positive by OnTrak TesTstik for Methamphetamine.

Clinical correlation for all positive and negative specimens is shown in the tables below.

TesTstik for Methamphetamine Assay Clinical Correlation (Cutoff = 500 ng/mL)

	Syva Emit II Negative	GC/MS positive (ng/mL)			Correlation	
		Near Cutoff		637- >10,000		
		396-497	525-614			
TesTstik for Methamphetamine	+	0	7	15	55	97%
	-	120*	8	0	0	

* Twelve percent of the screened negative clinical specimens were also confirmed negative by GC/MS.

Precision

The precision of OnTrak TesTstik for Methamphetamine was determined by testing 100 replicates per six various drug/drug metabolite concentration levels of multianalyte urine standards on three lots of OnTrak TesTstik for Methamphetamine (1800 tests in total). All lots are required to minimally perform with a greater than 95% confidence level that negative results will be attained with drugs at 25% ($\frac{1}{4}$ x) of the cutoff concentration, and greater than 95% confidence level that positive results will be attained with drugs at 150% ($1\frac{1}{2}$ x) of the cutoff concentration. The results are summarized in the table below.

X = Immunoassay Cutoff	Urine Standard Methamphetamine (ng/mL)	Precision Results (% Positive)
0 x	0	0
$\frac{1}{4}$ x	125	0
$\frac{1}{2}$ x	250	0
$\frac{3}{4}$ x	375	10
$1\frac{1}{4}$ x	625	100
$1\frac{1}{2}$ x	750	100

Specificity

The following structurally similar compounds for methamphetamine were tested for cross-reactivity using OnTrak TesTstik for Methamphetamine. Each compound tested was prepared in normal human urine. The results are expressed as that amount of the compound capable of giving a result equivalent to the cutoff.

Methamphetamine- related compounds	Approx. ng/mL equivalent to 500 ng/mL methamphetamine	Approx. percent cross-reactivity
<i>d</i> -Methamphetamine	500	100
<i>d,l</i> -Methamphetamine	600	83
3,4-Methylenedioxy- methamphetamine (MDMA)	1000	50
<i>l</i> -Methamphetamine	3500	14.3
Propylhexedrine	4000	12.5
<i>p</i> -Hydroxymethamphetamine	8000	6.25
Fenfluramine	10,000	5
3,4-Methylenedioxy-N- ethylamphetamine (MDEA)	20,000	2.5

1 Gently tear open the foil pouch to remove the stick. Pull the protective slide back to expose the sample pad.

2 Immerse the sample pad into the specimen until the specimen reaches the DIP LINE. Hold the stick in the specimen for approximately 10 seconds.

3 Remove the stick from the specimen. While holding the stick over the specimen cup, push the protective slide back over the sample pad.


Test complete—
ready to interpret
results

Test NOT complete—
DO NOT interpret
results

4 Allow the test to proceed until a blue band forms in the "TEST VALID" window (see example on the LEFT). Test is NOT complete until a blue band forms in the "TEST VALID" window (see example on the RIGHT). Once the "TEST VALID" line appears, remove the cover tab and interpret results.

5 A blue band in the result window indicates a NEGATIVE result.

6 No color in the result window indicates a POSITIVE result. The window appears as a white plus (+) sign.

<i>d,l</i> -Amphetamine	50,000	1
<i>d</i> -Amphetamine	50,000	1
Ranitidine	50,000	1
<i>l</i> -Phenylephrine	60,000	0.8
3,4-Methylenedioxy-amphetamine (MDA)	70,000	0.7
<i>l</i> -Amphetamine	100,000	0.5
<i>d,l</i> -Ephedrine	100,000	0.5
<i>l</i> -Ephedrine	100,000	0.5
β -Phenethylamine	100,000	0.5
<i>d</i> -Ephedrine	>100,000	<0.5
<i>p</i> -Hydroxyamphetamine	>100,000	<0.5
3-Hydroxytyramine	>100,000	<0.5
<i>l</i> -Norpseudoephedrine	>100,000	<0.5
Phentermine	>100,000	<0.5
<i>d,l</i> -Phenylpropanolamine	>100,000	<0.5
<i>d</i> -Phenylpropanolamine	>100,000	<0.5
<i>l</i> -Phenylpropanolamine	>100,000	<0.5
<i>d,l</i> -Pseudoephedrine	>100,000	<0.5
<i>l</i> -Pseudoephedrine	>100,000	<0.5
Tyramine	>100,000	<0.5
(4-Hydroxyphenethylamine)		

Cross-Reactivity with Additional Drugs

The following compounds were tested at 100,000 ng/mL and found not to cross-react with the OnTrak TesTstik for Methamphetamine assay.

Acetaminophen	Cocaine
Acetylsalicylic acid	Codeine
Aminopyrine	Cyclizine
Amitriptyline	Desipramine
Amobarbital	Dextromethorphan
Ampicillin	Dextropropoxyphene
Ascorbic acid	Diazepam
Aspartame	Diphenhydramine
Atropine	Diphenylhydantoin
Benzocaine	Dopamine
Benzoyllecgonine	Doxylamine
Benzphetamine	Ecgonine
Brompheniramine	Ecgonine methyl ester
Butabarbital	Epinephrine
Caffeine	Erythromycin
Calcium hypochlorite	Estril
Chlordiazepoxide	17 α -Ethinylestradiol
Chlorpheniramine	Fenoprofen
Chlorpromazine	Furosemide
Clemastine	Gentisic acid

Glutethimide	Oxazepam
Guaiacol glycerol ether	Penicillin G
Hydrochlorothiazide	Pentobarbital
Ibuprofen	Phenobarbital
Imipramine	Phenothiazine
Isoproterenol	Phenylbutazone
Ketamine	Phenyltoloxamine
LSD	Procyclidine
Lidocaine	Promethazine
Melanin	Quinidine
Meperidine	Quinine
Methadone	Secobarbital
Methapyrilene	Sulindac
Methaqualone	Tetracycline
Methylphenidate	Δ^9 -THC-9-carboxylic acid
Methypylon	Tetrahydrozoline
Morphine sulfate	Thioridazine
Naloxone	Trifluoperazine
Naltrexone	Trihexyphenidyl
Naproxen	Trimipramine
Niacinamide	Tripelennamine
Nordiazepam	Verapamil
Norethindrone	Zomepirac
Nortriptyline	

References

- Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
- Hoffman B, Lefkowitz R. Catecholamines, Sympathomimetic Drugs, and Adrenergic Receptor Antagonists. In: Hardman JG, Limbird LE, Molinoff PB, Ruddon RW, Goodman LS, Gilman A, eds. *Goodman and Gilman's The Pharmacological Basis of Therapeutics*. 9th ed. New York, NY: McGraw Hill Pub Co; 1996:199-248.
- Cook CE, Jeffcoat AR, Sadler BM, Hill JM, et al. Pharmacokinetics of Oral Methamphetamine and Effects of Repeated Daily Dosing in Humans. *Drug Metab Disposition*. 1993;20(6):856-862.
- Moore K. Amphetamines/Sympathomimetic Amines. In: Levine B, ed. *Principles of Forensic Toxicology*. American Association for Clinical Chemistry. 1999:265-285.
- Karch SB. *Drug Abuse Handbook*. CRC Press LLC 1998.
- Cook CE. Pyrolytic characteristics, pharmacokinetics, and bioavailability of smoked heroin, cocaine, phencyclidine and methamphetamine. *NIDA Research Monograph*. 1991;115:6-23.
- Data on file at Varian, Inc.

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Literatur

1. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
2. Hoffman B, Lefkowitz R. Catecholamines, Sympathomimetic Drugs, and Adrenergic Receptor Antagonists. In: Hardman JG, Limbird LE, Molinoff PB, Ruddon RW, Goodman LS, Gilman A, eds. *Goodman and Gilman's The Pharmacological Basis of Therapeutics*. 9th ed. New York, NY: McGraw Hill Pub Co; 1996:199-248.
3. Cook CE, Jeffcoat AR, Sadler BM, Hill JM, et al. Pharmacokinetics of Oral Methamphetamine and Effects of Repeated Daily Dosing in Humans. *Drug Metab Disposition*. 1993;20(6):856-862.
4. Moore K. Amphetamines/Sympathomimetic Amines. In: Levine B, ed. *Principles of Forensic Toxicology*. American Association for Clinical Chemistry. 1999:265-285.
5. Karch SB. *Drug Abuse Handbook*. CRC Press LLC 1998.
6. Cook CE. Pyrolytic characteristics, pharmacokinetics, and bioavailability of smoked heroin, cocaine, phencyclidine and methamphetamine. *NIDA Research Monograph*. 1991;115:6-23.
7. Dokumentation, Varian, Inc.

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Bibliographie

- Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
- Hoffman B, Lefkowitz R. Catecholamines, Sympathomimetic Drugs, and Adrenergic Receptor Antagonists. In: Hardman JG, Limbird LE, Molinoff PB, Ruddon RW, Goodman LS, Gilman A, eds. *Goodman and Gilman's The Pharmacological Basis of Therapeutics*. 9th ed. New York, NY: McGraw Hill Pub Co; 1996:199-248.
- Cook CE, Jeffcoat AR, Sadler BM, Hill JM, et al. Pharmacokinetics of Oral Methamphetamine and Effects of Repeated Daily Dosing in Humans. *Drug Metab Disposition*. 1993;20(6):856-862.
- Moore K. Amphetamines/Sympathomimetic Amines. In: Levine B, ed. *Principles of Forensic Toxicology*. American Association for Clinical Chemistry. 1999:265-285.
- Karch SB. *Drug Abuse Handbook*. CRC Press LLC 1998.
- Cook CE. Pyrolytic characteristics, pharmacokinetics, and bioavailability of smoked heroin, cocaine, phencyclidine and methamphetamine. *NIDA Research Monograph*. 1991;115:6-23.
- Documentation de Varian, Inc.

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- Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
- Hoffman B, Lefkowitz R. Catecholamines, Sympathomimetic Drugs, and Adrenergic Receptor Antagonists. En: Hardman JG, Limbird LE, Molinoff PB, Ruddon RW, Goodman LS, Gilman A, eds. *Goodman and Gilman's The Pharmacological Basis of Therapeutics*. 9th ed. New York, NY: McGraw Hill Pub Co; 1996:199-248.
- Cook CE, Jeffcoat AR, Sadler BM, Hill JM, et al. Pharmacokinetics of Oral Methamphetamine and Effects of Repeated Daily Dosing in Humans. *Drug Metab Disposition*. 1993;20(6):856-862.
- Moore K. Amphetamines/Sympathomimetic Amines. En: Levine B, ed. *Principles of Forensic Toxicology*. American Association for Clinical Chemistry. 1999:265-285.
- Karch SB. *Drug Abuse Handbook*. CRC Press LLC 1998.
- Cook CE. Pyrolytic characteristics, pharmacokinetics, and bioavailability of smoked heroin, cocaine, phencyclidine and methamphetamine. *NIDA Research Monograph*. 1991;115:6-23.
- Datos de archivo, Varian, Inc.

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Riferimenti

1. Hawks RL, Chiang CN (a cura di) Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
2. Hoffman B, Lefkowitz R. Catecholamines, Sympathomimetic Drugs, and Adrenergic Receptor Antagonists. In: Hardman JG, Limbird LE, Molinoff PB, Ruddon RW, Goodman Gilman A (a cura di) *Goodman and Gilman's The Pharmacological Basis of Therapeutics*. nona ed. New York, NY: McGraw Hill Pub Co; 1996:199-248.
3. Cook CE, Jeffcoat AR, Sadler BM, Hill JM, et al. Pharmacokinetics of Oral Methamphetamine and Effects of Repeated Daily Dosing in Humans. *Drug Metab Disposition*. 1993;20(6):856-862.
4. Moore K. Amphetamines/Sympathomimetic Amines. In: Levine B (a cura di) *Principles of Forensic Toxicology*. American Association for Clinical Chemistry. 1999:265-285.
5. Karch SB. *Drug Abuse Handbook*. CRC Press LLC 1998.
6. Cook CE. Pyrolytic characteristics, pharmacokinetics, and bioavailability of smoked heroin, cocaine, phencyclidine and methamphetamine. *NIDA Research Monograph*. 1991;115:6-23.
7. Dati presenti nel file, Varian, Inc.

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Bibliografia

- Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
- Hoffman B, Lefkowitz R. Catecholamines, Sympathomimetic Drugs, and Adrenergic Receptor Antagonists. En: Hardman JG, Limbird LE, Molinoff PB, Ruddon RW, Goodman Gilman A, eds. *Goodman and Gilman's The Pharmacological Basis of Therapeutics*. 9th ed. New York, NY: McGraw Hill Pub Co; 1996:199-248.
- Cook CE, Jeffcoat AR, Sadler BM, Hill JM, et al. Pharmacokinetics of Oral Methamphetamine and Effects of Repeated Daily Dosing in Humans. *Drug Metab Disposition*. 1993;20(6):856-862.
- Moore K. Amphetamines/Sympathomimetic Amines. En: Levine B, ed. *Principles of Forensic Toxicology*. American Association for Clinical Chemistry. 1999:265-285.
- Karch SB. *Drug Abuse Handbook*. CRC Press LLC 1998.
- Cook CE. Pyrolytic characteristics, pharmacokinetics, and bioavailability of smoked heroin, cocaine, phencyclidine and methamphetamine. *NIDA Research Monograph*. 1991;115:6-23.
- Ficheiros da Varian, Inc.

www.varian-onsite.com

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Dentro dos E.U.A.: 1-800-737-9667

Fora dos E.U.A.: 00 1 (949) 770-9381, extensão 2

DAT.TechSupport@varianinc.com

 Varian, Inc.
Lake Forest, CA 92630 USA

 Varian B.V.
Herculesweg 8, 4338 PL Middelburg
 Holanda

Número da patente: 5,770,458

16141 1

© 2003, Varian, Inc.

VARIAN