English Intended Use

OnTrak TesTstik™ for Morphine is an *in vitro* diagnostic test intended for professional use for the qualitative detection of morphine in urine at or above a cutoff concentration of 300 ng/mL. OnTrak TesTstik for Morphine is not intended for over-the-counter sale.

OnTrak TesTstik provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmation method.¹ Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Summary

Morphine, a natural product of the opium poppy, is a narcotic analgesic used for centuries as a medicine for the relief of severe pain. Extracted from opium obtained from the poppy's resin, morphine may, in turn, be further chemically refined to heroin (the more potent, diacetylated analog of the parent drug). These chemically similar "opiates" reduce sensitivity to physical and psychological stimuli, dulling pain, fear and anxiety. Users are usually lethargic and indifferent. Accompanying effects may include constriction of the pupils, itching, constipation, nausea, vomiting, and respiratory depression. Death by overdose, usually resulting from dose miscalculation or dose-strength variability, is caused by respiratory failure.²⁻⁴

The opiates are usually administered intravenously or subcutaneously, but may also be smoked or sniffed. Upon entering the circulation, they tend to concentrate in the lungs, spleen, kidneys, and liver; lower concentrations are found in the body's musculature and central nervous system. A variety of pathways are involved in the body's detoxification of the opiates, including the removal of chemical side groups (dealkylation), addition of hydroxyl groups, hydrolytic breakdown, and conjugation to glucuronic acid (a common body sugar).⁵ Morphine is excreted in the urine as morphine-3-glucuronide, unchanged free morphine, and other minor metabolites. Although some opiate metabolites appear in the bile and feces, urinary excretion is the primary route of elimination.^{1,6}

The opiates produce strong physical dependence; withdrawal symptoms can begin to appear within a few hours of the last dose and may continue for 5 - 10 days. The addict may pursue continued opiate use as much to avoid the discomfort of withdrawal as to achieve the desired insensate euphoria.^{7,8}

Test principle

OnTrak TesTstik assays are based on the principles of microparticle capture inhibition. The test relies on the competition between drug, which may be present in the urine being tested, and drug conjugate immobilized on membrane for binding to antibody coated onto colored microparticles. When the TesTstik is immersed in the urine specimen, some of the specimen is absorbed into the TesTstik sample pad. The absorbed sample travels through a reagent strip contained in the device by capillary action. In the reagent strip, the specimen rehydrates and mobilizes antibody-coated blue microparticles. The microparticle-urine suspension continues to migrate through the reagent strip and comes in contact with the immobilized drug conjugate. In the absence of drug in the urine, the antibody-coated microparticles bind to the drug conjugate and a blue band is formed at the result window ("negative" sign).

When drug is present in the specimen, it binds to the antibody-coated microparticles. If sufficient drug is present, the microparticles are inhibited from binding the drug conjugate and no blue band is formed at the result window. A presumptive positive ("non-negative") sample causes the membrane to remain white ("positive" sign).

An additional antibody/antigen reaction occurs at the "TEST VALID" area. The "TEST VALID" blue band forms when anti-bovine serum albumin (anti-BSA) antibodies, which are imbedded in the reagent membrane, bind to the BSA antigen on the blue microparticles. The presence of the "TEST VALID" band indicates that the test has completed, the reagents are viable, and, the results are ready to interpret.

Reagents

1. Blue-dyed microparticles coated with mouse monoclonal anti-morphine antibody, and BSA in a buffered solution containing preservative and dried onto a membrane.
2. Drug conjugates immobilized on membrane.
3. Mouse monoclonal anti-BSA antibody immobilized on membrane.

Precautions and Warnings

For *in vitro* diagnostic use. Exercise the normal precautions for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines.

Storage and Stability

OnTrak TesTstik should be stored at room temperature, 15 - 30 °C (59 - 86 °F), in the original foil pouch. If the foil pouch containing the TesTstik is damaged (e.g., a hole or tear), do not use that TesTstik. Do not use the OnTrak TesTstik after the expiration date printed on the foil pouch.

Specimen Collection and Preparation

OnTrak TesTstik is formulated for use with urine specimens. Fresh urine specimens do not require any special handling or pretreatment. No additives or preservatives are required. It is recommended that specimens be stored refrigerated (2 - 8 °C or 36 - 46 °F) and tested within two days of collection. For prolonged storage, freezing is recommended.

Specimens containing human-sourced materials should be handled as if potentially infectious. Use safe laboratory procedures such as those outlined in *Biosafety in Microbiological and Biomedical Laboratories* (HHS Publication Number [CDC] 93-8395).

Testing Procedure*Materials provided*

50 individually packaged OnTrak TesTstik for Morphine sticks

Ancillary materials

	Part No.	Quantity
Evidence Tape	42212	100 strips/pkg.
Latex Gloves	42210	50 pairs
Collection Cups	47123	100 cups/box
Adulteration Strips	3111407	25 strips/box

Contact Customer Technical Support for information regarding commercially available control material.

Assay

OnTrak TesTstik is ready for use directly from the foil pouch and must be used within eight hours after removal from the pouch. This time frame may be shorter in high humidity climates. Do not open the foil pouch unless it is at room temperature (15 - 30 °C or 59 - 86 °F).

Specimens should be between room temperature and body temperature 15 - 37 °C (59 - 99 °F). If the specimen was refrigerated, allow it to reach room temperature before testing. To ensure a sufficient quantity of sample for confirmation testing (as necessary), collect a minimum of 30 mL of sample.

1. To expose the sample pad, pull the protective slide back until it stops (up to the "Dip line").
2. Gently immerse the exposed sample pad into the specimen until the specimen reaches the dip line. Do not immerse the entire device into the specimen. Keep the OnTrak TesTstik for Morphine in the specimen for approximately 10 seconds.
3. Remove the OnTrak TesTstik for Morphine from the specimen and while holding the OnTrak TesTstik for Morphine over the specimen, push the protective slide back over the sample pad until it stops. Allow the test to proceed until a blue band is observed in the "TEST VALID" window (see Results section below). The OnTrak TesTstik for Morphine may be placed down or held while the test is developing. There is no required timing step. Note: If the Test Valid band does not appear, re-dip the OnTrak TesTstik for Morphine for approximately 10 seconds. Allow to sit undisturbed for an additional 1-3 minutes. Once the Test Valid band forms, proceed to step 4.
4. Remove the tab that covers the result window from the protective slide. Read and record the result.

Quality Control

Quality control testing at regular intervals is good laboratory practice and may be required by federal, state, or local guidelines. Always check with the appropriate licensing or accrediting bodies to ensure your quality control program meets the established standards.

Internal control: Each OnTrak TesTstik MOR has built-in process controls. The "TEST VALID" band is an internal POSITIVE process control. A blue "TEST VALID" band should always appear if an adequate sample volume is used, the sample and reagent are wicking on the membrane, and the reagents in the "TEST VALID" band and the conjugate-color indicator are working. In addition, the background in the result window should become clear and provide a distinct result. If smears of color appear in the background, nonspecific binding of the antibody may be occurring, making the test result difficult to interpret. The clearing of the result window background may be considered an internal NEGATIVE process control. If the "TEST VALID" band does not appear in the test window, or the background does not clear, then the test is not valid and a new test must be performed.

External control: The use of external positive and negative controls is recommended to test each shipment of product or when a new lot is used. Contact Customer Technical Support for information regarding commercially available control material. The controls should produce the expected results. If expected results are not obtained, do not run test specimens, and contact Customer Technical Support. Follow the appropriate federal, state and local guidelines when running external controls.

Results

Once a full, partial or broken blue band is formed in the "TEST VALID" window (in approximately 3 minutes) remove the cover tab and interpret results as either presumptive positive or negative. Results may be interpreted for up to 30 minutes after the "TEST VALID" band forms.

A negative result is the presence of a blue band ("negative" sign) anywhere in the test result window. Since a negative sample may also give a faint or incomplete color band, **any** full, partial or broken band in the test result window indicates a negative result.

A presumptive positive result (drug present at or above the cutoff) is the lack of a blue band in the test result window. The window appears as a white plus (+) sign.

Limitations

See Specific Performance Characteristics for information on substances tested for cross-reactivity in this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors). Adulteration of the sample with oxidizing agents and/or acidic or alkaline agents can cause erroneous results. If adulteration is suspected, another sample

should be collected and a new OnTrak TesTstik used. For information on how adulterants may affect the performance of this assay, contact Customer Technical Support.

A presumptive positive result with this assay indicates the presence of drug or drug metabolites in urine but does not reflect the degree of intoxication.

Ingestion of food containing poppy seeds may cause positive morphine test results due to the natural content of morphine in poppy seeds.⁹ When present in sufficient quantity, this "poppy seed morphine" can be detected by all chemical and immunological assays for morphine.

Specific Performance Characteristics¹⁰

Accuracy

OnTrak TesTstik for Morphine was evaluated using specimens screened by an automated immunoassay and confirmed positive by GC/MS (at a 300 ng/mL cutoff). Forty nine samples (>99%) positive for morphine were positive by OnTrak TesTstik.

		GC/MS	
		+	-
TesTstik	+	49	
	-	1*	

*Sample contained 394 ng/mL morphine by GC/MS.

One hundred six urine samples, obtained from a clinical laboratory and screened negative by an automated immunoassay relative to a 300 ng/mL cutoff for morphine were evaluated using OnTrak TesTstik. One hundred five (>99%) were negative for morphine by OnTrak TesTstik. All positive and negative samples were also assayed by, and compared to, Abuscreen OnTrak for Morphine. One hundred fifty six samples tested by both OnTrak TesTstik and Abuscreen OnTrak for Morphine demonstrated 99.3% agreement.

		Abuscreen OnTrak	
		+	-
TesTstik	+	50	0
	-	1*	105

*Sample contained 394 ng/mL morphine by GC/MS.

Precision

The precision of OnTrak TesTstik was determined by testing various concentrations of drug standards for the primary drugs or drug metabolites detected by this system. All lots are required to minimally perform with a greater than 95% confidence level that negative results will be attained with drugs at 25% ($\frac{1}{4}$ x) of their respective cutoff concentrations, and a greater than 95% confidence level that positive results will be attained with drugs at 150% ($1\frac{1}{2}$ x) of their respective cutoff concentrations.

Specificity

The following structurally similar compounds for morphine were tested for cross-reactivity with OnTrak TesTstik for Morphine. Each compound tested was prepared in normal human urine. The results are expressed as that amount of the compound capable of giving a result equivalent to the cutoff.

Morphine-related compounds	Approx. ng/mL equivalent to 300 ng/mL morphine	Approx. percent cross-reactivity
Codeine	300	100
Ethyl morphine HCl	300	100
Dihydrocodeine bitartrate	300	100
Morphine-3-glucuronide	400	75
Dihydromorphine	400	75
6-Acetyl-morphine	400	75
Hydrocodone bitartrate	600	50
Hydromorphone HCl	700	43
Thebaine	2000	15
N-Norcodeine HCl	25,000	1.2
Oxycodone	25,000	1.2
Meperidine	100,000	0.3

1 Gently tear open the foil pouch to remove the stick. Pull the protective slide back to expose the sample pad.

2 Immerse the sample pad into the specimen until the specimen reaches the DIP LINE. Hold the stick in the specimen for approximately 10 seconds.

3 Remove the stick from the specimen. While holding the stick over the specimen cup, push the protective slide back over the sample pad.


Test complete— ready to interpret results

Test NOT complete— DO NOT interpret results

4 Allow the test to proceed until a blue band forms in the “TEST VALID” window (see example on the LEFT). Test is NOT complete until a blue band forms in the “TEST VALID” window (see example on the RIGHT). Once the “TEST VALID” line appears, remove the cover tab and interpret results.

5 A blue band in the result window indicates a NEGATIVE result.

6 No color in the result window indicates a POSITIVE result. The window appears as a white plus (+) sign.

Cross-Reactivity with Additional Drugs

The following compounds were tested at 100,000 ng/mL and found not to cross-react with the OnTrak TestStik for Morphine assay.

Acetaminophen	Fenopropfen
Acetylsalicylic acid	Furosemide
Aminopyrine	Gentisic acid
Amitriptyline	Glutethimide
Amobarbital	Guaiacol glycerol ether
Amphetamine	Hydrochlorothiazide
Ampicillin	p-Hydroxyamphetamine
Ascorbic acid	Ibuprofen
Aspartame	Imipramine
Atropine	Isoproterenol
Benzocaine	Ketamine
Benzoylcegonine	LSD
Benzphetamine	Lidocaine
Butabarbital	MDA
Caffeine	MDMA
Calcium hypochlorite	Melanin
Chlordiazepoxide	Methadone
Chloroquine	Methamphetamine
Chlorpheniramine	Methaqualone
Chlorpromazine	Methylphenidate
Cocaine	Methyprylon
Dextromethorphan	Naloxone
Dextropropoxyphene	Naltrexone
Diazepam	Naproxen
Diphenhydramine	Niacinamide
Diphenylhydantoin	Nordiazepam
Dopamine	Norethindrone
Doxylamine	Norpseudoephedrine
Ecgonine	Oxazepam
Ecgonine methyl ester	Penicillin G
Ephedrine	Pentobarbital
Epinephrine	Phencyclidine
Erythromycin	β-Phenethylamine
Estriol	Phenobarbital

Phenothiazine	Quinine
Phentermine	Secobarbital
Phenylbutazone	Sulindac
Phenylephrine	Tetracycline
Phenylpropranolamine	Δ ⁹ -THC-9-carboxylic acid
Procaine	Tetrahydrozoline
Promethazine	Trifluoperazine
Pseudoephedrine	Tyramine
Quinidine	Verapamil

References

1. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
2. Papadopoulos CN, Keats AS. Studies of analgesic drugs VI: Comparative respiratory depressant activity of phenazocine and morphine. *Clin Pharmacol Ther*. 1961;2:8.
3. Louria DB, et al. The major medical complications of heroin addiction. *Ann Intern Med*. 1967;67:1.
4. Dripps RD, Comroe JH Jr. Clinical studies on morphine I: The immediate effect of morphine administered intravenously and intramuscularly upon the respiration of normal man. *Anesthesiology*. 1945;6:462.
5. Way EL, Adler TK. The pharmacological implications of the fate of morphine and its surrogates. *Pharmacol Rev*. 1960;12:383.
6. Boerner U, et al. The metabolism of morphine and heroin in man. *Metab Rev*. 1975;4:39.
7. Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc; 1984.
8. Jaffe JH. Drug addiction and drug abuse. in: Gilman A, Rall T, Neis a, Taylor P, eds. *Goodman and Gillman's The Pharmacological Basis of Therapeutics*. 8th ed. New York, NY: McGraw Hill Pub Co; 1993:187-220.
9. EISohly HN, EISohly MA. Poppy seed ingestion and opiates urinalysis: A closer look. *J Anal Toxicol*. 1990;14:308-310.
10. Data on file at Varian, Inc.

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Literatur

- Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
- Papadopoulos CN, Keats AS. Studies of analgesic drugs VI: Comparative respiratory depressant activity of phenazocine and morphine. *Clin Pharmacol Ther.* 1961;2:8.
- Louria DB, et al. The major medical complications of heroin addiction. *Ann Intern Med.* 1967;67:1.
- Dripps RD, Comroe JH Jr. Clinical studies on morphine I: The immediate effect of morphine administered intravenously and intramuscularly upon the respiration of normal man. *Anesthesiology.* 1945;6:462.
- Way EL, Adler TK. The pharmacological implications of the fate of morphine and its surrogates. *Pharmacol Rev.* 1960;12:383.
- Boerner U, et al. The metabolism of morphine and heroin in man. *Metab Rev.* 1975;4:39.
- Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc; 1984.
- Jaffe JH. Drug addiction and drug abuse. in: Gilman A, Rall T, Neis a, Taylor P, eds. *Goodman and Gillman's The Pharmacological Basis of Therapeutics*. 8th ed. New York, NY: McGraw Hill Pub Co; 1993:187-220.
- EISohly HN, EISohly MA. Poppy seed ingestion and opiates urinalysis: A closer look. *J Anal Toxicol.* 1990;14:308-310.
- Dokumentation, Varian, Inc.

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Bibliographie

- Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
- Papadopoulos CN, Keats AS. Studies of analgesic drugs VI: Comparative respiratory depressant activity of phenazocine and morphine. *Clin Pharmacol Ther.* 1961;2:8.
- Louria DB, et al. The major medical complications of heroin addiction. *Ann Intern Med.* 1967;67:1.
- Dripps RD, Comroe JH Jr. Clinical studies on morphine I: The immediate effect of morphine administered intravenously and intramuscularly upon the respiration of normal man. *Anesthesiology.* 1945;6:462.
- Way EL, Adler TK. The pharmacological implications of the fate of morphine and its surrogates. *Pharmacol Rev.* 1960;12:383.
- Boerner U, et al. The metabolism of morphine and heroin in man. *Metab Rev.* 1975;4:39.
- Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc; 1984.
- Jaffe JH. Drug addiction and drug abuse. in: Gilman A, Rall T, Neis a, Taylor P, eds. *Goodman and Gillman's The Pharmacological Basis of Therapeutics*. 8th ed. New York, NY: McGraw Hill Pub Co; 1993:187-220.
- EISohly HN, EISohly MA. Poppy seed ingestion and opiates urinalysis: A closer look. *J Anal Toxicol.* 1990;14:308-310.
- Documentation de Varian, Inc.

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- Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
- Papadopoulos CN, Keats AS. Studies of analgesic drugs VI: Comparative respiratory depressant activity of phenazocine and morphine. *Clin Pharmacol Ther.* 1961;2:8.
- Louria DB, et al. The major medical complications of heroin addiction. *Ann Intern Med.* 1967;67:1.
- Dripps RD, Comroe JH Jr. Clinical studies on morphine I: The immediate effect of morphine administered intravenously and intramuscularly upon the respiration of normal man. *Anesthesiology.* 1945;6:462.
- Way EL, Adler TK. The pharmacological implications of the fate of morphine and its surrogates. *Pharmacol Rev.* 1960;12:383.
- Boerner U, et al. The metabolism of morphine and heroin in man. *Metab Rev.* 1975;4:39.
- Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc; 1984.
- Jaffe JH. Drug addiction and drug abuse. En: Gilman A, Rall T, Neis a, Taylor P, eds. *Goodman and Gillman's The Pharmacological Basis of Therapeutics*. 8th ed. New York, NY: McGraw Hill Pub Co; 1993:187-220.
- EISohly HN, EISohly MA. Poppy seed ingestion and opiates urinalysis: A closer look. *J Anal Toxicol.* 1990;14:308-310.
- Datos de archivo, Varian, Inc.

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Riferimenti

- Hawks RL, Chiang CN (a cura di) Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
- Papadopoulos CN, Keats AS. Studies of analgesic drugs VI: Comparative respiratory depressant activity of phenazocine and morphine. *Clin Pharmacol Ther.* 1961;2:8.
- Louria DB, et al. The major medical complications of heroin addiction. *Ann Intern Med.* 1967;67:1.
- Dripps RD, Comroe JH Jr. Clinical studies on morphine I: The immediate effect of morphine administered intravenously and intramuscularly upon the respiration of normal man. *Anesthesiology.* 1945;6:462.
- Way EL, Adler TK. The pharmacological implications of the fate of morphine and its surrogates. *Pharmacol Rev.* 1960;12:383.
- Boerner U, et al. The metabolism of morphine and heroin in man. *Metab Rev.* 1975;4:39.
- Blum K. *Handbook of Abusable Drugs*. prima ed. New York, NY: Gardner Press Inc; 1984.
- Jaffe JH. Drug addiction and drug abuse. in: Gilman A, Rall T, Neis A, Taylor P (a cura di) *Goodman and Gillman's The Pharmacological Basis of Therapeutics*. ottava ed. New York, NY: McGraw Hill Pub Co; 1993:187-220.
- El Sohly HN, El Sohly MA. Poppy seed ingestion and opiates urinalysis: A closer look. *J Anal Toxicol.* 1990;14:308-310.
- Dati presenti nel file, Varian, Inc.

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Bibliografia

- Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
- Papadopoulos CN, Keats AS. Studies of analgesic drugs VI: Comparative respiratory depressant activity of phenazocine and morphine. *Clin Pharmacol Ther.* 1961;2:8.
- Louria DB, et al. The major medical complications of heroin addiction. *Ann Intern Med.* 1967;67:1.
- Dripps RD, Comroe JH Jr. Clinical studies on morphine I: The immediate effect of morphine administered intravenously and intramuscularly upon the respiration of normal man. *Anesthesiology.* 1945;6:462.
- Way EL, Adler TK. The pharmacological implications of the fate of morphine and its surrogates. *Pharmacol Rev.* 1960;12:383.
- Boerner U, et al. The metabolism of morphine and heroin in man. *Metab Rev.* 1975;4:39.
- Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc; 1984.
- Jaffe JH. Drug addiction and drug abuse. in: Gilman A, Rall T, Neis a, Taylor P, eds. *Goodman and Gilman's The Pharmacological Basis of Therapeutics*. 8th ed. New York, NY: McGraw Hill Pub Co; 1993:187-220.
- EISohly HN, EISohly MA. Poppy seed ingestion and opiates urinalysis: A closer look. *J Anal Toxicol.* 1990;14:308-310.
- Ficheiros da Varian, Inc.

www.varian-onsite.com

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Dentro dos E.U.A.: 1-800-737-9667

Fora dos E.U.A.: 00 1 (949) 770-9381, extensão 2

DAT.TechSupport@varianinc.com

 Varian, Inc.
Lake Forest, CA 92630 USA

 Varian B.V.
Herculesweg 8, 4338 PL Middelburg
 Holanda

Número da patente: 5,770,458

16152 1

© 2003, Varian, Inc.

VARIAN