English

## Intended Use

OnTrak TesTstik™ for Phencyclidine (PCP) is an *in vitro* diagnostic test intended for professional use for the qualitative detection of PCP in urine. The OnTrak TesTstik for PCP cutoff level (25 ng/mL) is based on US federal mandatory guidelines.<sup>1</sup> OnTrak TesTstik for Phencyclidine (PCP) is not intended for over-the-counter sale.

**OnTrak TesTstik provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmation method.<sup>2</sup> Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.**

## Summary

Phencyclidine (PCP) or "Angel Dust" is an arylcyclohexylamine that is structurally related to the anesthetic ketamine and was originally developed as an intravenous analgesic/anesthetic intended for human use. However, the occurrence of side effects symptomatic of emergence psychosis negated its potential clinical utility. It was then made available as a veterinary tranquilizer/anesthetic.<sup>3-5</sup>

PCP also possesses hallucinogenic, central nervous system (CNS)-stimulant, CNS-depressant, and analgesic properties, the expression of which is dose- and species-dependent.<sup>5</sup> Because of its simple and inexpensive chemical synthesis, PCP is available in many geographic areas as a substitute for mescaline and LSD.<sup>3</sup> The water-soluble powder can be ingested, inhaled, or injected intravenously, but it is most commonly smoked.<sup>6</sup> PCP is also commonly used in combination with other illicit drugs, as are its structurally related analogs: TCP, PHP, PCE, and ketamine.<sup>3,5</sup> Typical street doses (1 - 10 mg) can cause tachycardia, hypertension, hallucinations, stupor, disorientation, lethargy, sensory isolation, and gross loss of coordination. Excitation and agitation may also occur, leading to unpredictably violent behavior not usually encountered with other hallucinogens. Higher doses can culminate in convulsions and prolonged or fatal coma. Neither the development of tolerance to PCP nor the existence of a withdrawal syndrome have been firmly established. Physical dependence on this drug, similar to that seen with the opiates and other CNS depressants, has not been observed.<sup>3</sup>

PCP is readily absorbed across mucous membranes into the blood, and its immediate effects are relatively short-lived. It is metabolized by the liver via ring hydroxylation to two hydroxylated metabolites and a dihydroxylated species that are excreted as glucuronide conjugates in the urine.<sup>7,8</sup>

## Test Principle

OnTrak TesTstik assays are based on the principles of microparticle capture inhibition. The test relies on the competition between drug, which may be present in the urine being tested, and drug conjugate immobilized on membrane for binding to antibody coated onto colored microparticles. When the TesTstik is immersed in the urine specimen, some of the specimen is absorbed into the TesTstik sample pad. The absorbed specimen travels through a reagent strip contained in the device by capillary action. In the reagent strip, the specimen rehydrates and mobilizes antibody-coated blue microparticles. The microparticle-urine suspension continues to migrate through the reagent strip and comes in contact with the immobilized drug conjugate. In the absence of drug in the urine, the antibody-coated microparticles bind to the drug conjugate and a blue band is formed at the result window ("negative" sign).

When drug is present in the specimen, it binds to the antibody-coated microparticles. If sufficient drug is present, the microparticles are inhibited

from binding the drug conjugate and no blue band is formed at the result window. A presumptive positive ("non-negative") sample causes the membrane to remain white ("positive" sign).

An additional antibody/antigen reaction occurs at the "TEST VALID" area. The "TEST VALID" blue band forms when anti-bovine serum albumin (anti-BSA) antibodies, which are imbedded in the reagent membrane, bind to the BSA antigen on the blue microparticles. The presence of the "TEST VALID" band indicates that the test has completed, the reagents are viable, and the results are ready to interpret.

## Reagents

1. Blue-dyed microparticles coated with mouse monoclonal anti-phencyclidine antibody, and BSA in a buffered solution containing preservative and dried onto a membrane.
2. Drug conjugates immobilized on membrane.
3. Mouse monoclonal anti-BSA antibody immobilized on membrane.

## Precautions and Warnings

For *in vitro* diagnostic use. Exercise the normal precautions for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines.

## Storage and stability

OnTrak TesTstik should be stored at room temperature, 59 - 86 °F (15 - 30 °C), in the original foil pouch. If the foil pouch containing the TesTstik is damaged (e.g., a hole or tear), do not use that TesTstik. Do not use the TesTstik after the expiration date printed on the foil pouch.

## Specimen Collection and Preparation

OnTrak TesTstik is formulated for use with urine specimens. Fresh urine specimens do not require any special handling or pretreatment. No additives or preservatives are required. It is recommended that specimens be stored refrigerated (2 - 8 °C or 36 - 46 °F) and tested within two days of collection. For prolonged storage, freezing is recommended.

Specimens containing human-sourced materials should be handled as if potentially infectious. Use safe laboratory procedures such as those outlined in *Biosafety in Microbiological and Biomedical Laboratories* (HHS Publication Number [CDC] 93-8395).

## Testing Procedure

### Materials provided

50 individually packaged OnTrak TesTstik for PCP sticks

### Ancillary materials

	Part No.	Quantity
Evidence Tape	42212	100 strips/pkg.
Latex Gloves	42210	50 pairs
Collection Cups	47123	100 cups/box
Adulteration Strips	3111407	25 strips/box

Contact Customer Technical Support for information regarding commercially available control material.

## Assay

OnTrak TesTstik is ready for use directly from the foil pouch and must be used within eight hours after removal from the pouch. This time frame may be shorter in high humidity climates. Do not open the foil pouch unless it is at room temperature (15 - 30 °C or 59 - 86 °F).

Specimens should be between room temperature and body temperature 15 - 37 °C (59 - 99 °F). If the specimen was refrigerated, allow it to reach room temperature before testing. To ensure a sufficient quantity of sample for confirmation testing (as necessary), collect a minimum of 30 mL of sample.

- To expose the sample pad, pull the protective slide back until it stops (up to the "Dip" line).
- Gently immerse the exposed sample pad into the specimen until the specimen reaches the dip line. Do not immerse the entire device into the specimen. Keep the OnTrak TesTstik for PCP in the specimen for approximately 10 seconds.
- Remove the OnTrak TesTstik for PCP from the specimen and while holding the OnTrak TesTstik for PCP over the specimen, push the protective slide back over the sample pad until it stops. Allow the test to proceed until a blue band is observed in the "TEST VALID" window (see Results section below). The OnTrak TesTstik for PCP may be placed down or held while the test is developing. There is no required timing step. Note: If the Test Valid band does not appear, re-dip the OnTrak TesTstik for PCP for approximately 10 seconds. Allow to sit undisturbed for an additional 1-3 minutes. Once the Test Valid band forms, proceed to step 4.
- Remove the tab that covers the result window from the protective slide. Read and record the result.

### Quality Control

Quality control testing at regular intervals is good laboratory practice and may be required by federal, state, or local guidelines. Always check with the appropriate licensing or accrediting bodies to ensure your quality control program meets the established standards.

**Internal control:** Each OnTrak TesTstik PCP has built-in process controls. The "TEST VALID" band is an internal POSITIVE process control. A blue "TEST VALID" band should always appear if an adequate sample volume is used, the sample and reagent are wicking on the membrane, and the reagents in the "TEST VALID" band and the conjugate-color indicator are working. In addition, the background in the result window should become clear and provide a distinct result. If smears of color appear in the background, nonspecific binding of the antibody may be occurring, making the test result difficult to interpret. The clearing of the result window background may be considered an internal NEGATIVE process control. If the "TEST VALID" band does not appear in the test window, or the background does not clear, then the test is not valid and a new test must be performed.

**External control:** The use of external positive and negative controls is recommended to test each shipment of product or when a new lot is used. Contact Customer Technical Support for information regarding commercially available control material. The controls should produce the expected results. If expected results are not obtained, do not run test specimens, and contact Customer Technical Support. Follow the appropriate federal, state and local guidelines when running external controls.

### Results

Once a full, partial or broken blue band is formed in the "TEST VALID" window (in approximately 3 minutes) remove the cover tab and interpret results as either presumptive positive or negative. Results may be interpreted for up to 30 minutes after the "TEST VALID" band forms.

A negative result is the presence of a blue band ("negative" sign) anywhere in the test result window. Since a negative sample may give a faint or incomplete color band, **any** full, partial or broken band in the test result window indicates a negative result.

A presumptive positive result (drug present at or above the cutoff) is the lack of a blue band in the test result window. The window appears as a white plus (+) sign.

### Limitations

See Specific Performance Characteristics for information on substances tested for cross-reactivity in this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors).

Adulteration of the sample with oxidizing agents and/or acidic or alkaline agents can cause erroneous results. If adulteration is suspected, another sample should be collected and a new OnTrak TesTstik used. For information on how adulterants may affect the performance of this assay, contact Customer Technical Support.

A presumptive positive result with this assay indicates the presence of drug or drug metabolites in urine but does not reflect the degree of intoxication.

## Specific Performance Characteristics<sup>9</sup>

### Accuracy

OnTrak TesTstik for PCP was evaluated using specimens screened by an automated immunoassay and confirmed positive by GC/MS (at a 25 ng/mL cutoff). Fifty samples positive for PCP were positive by OnTrak TesTstik for PCP (100%).

		GC/MS	
		+	
TesTstik	+	50	
	-	0	

One hundred six urine samples, obtained from a clinical laboratory and screened negative by an automated immunoassay relative to a 25 ng/mL cutoff for PCP, were evaluated and found negative using OnTrak TesTstik for PCP.

All positive and negative samples were also assayed by, and compared to, Abuscreen ONLINE for PCP. All samples demonstrated 100% agreement between the two assays.

		Abuscreen OnTrak	
		+	-
TesTstik	+	50	0
	-	0	106

### Precision

The precision of OnTrak TesTstik for PCP was determined by testing various concentrations of drug standards for the primary drugs or drug metabolites detected by this system. All lots are required to minimally perform with a greater than 95% confidence level that negative results will be attained with drugs at 25% ( $\frac{1}{4}x$ ) of their respective cutoff concentrations, and a greater than 95% confidence level that positive results will be attained with drugs at 150% ( $\frac{1}{2}x$ ) of their respective cutoff concentrations.

### Specificity

The following structurally similar compounds for PCP were tested for cross-reactivity with OnTrak TesTstik for PCP. Each compound tested was prepared in normal human urine. The results are expressed as that amount of the compound capable of giving a result equivalent to the cutoff.

PCP-related compounds	Approx. ng/mL equivalent to 25 ng/mL phencyclidine	Approx. percent cross-reactivity
1-(1-Phenylcyclohexyl) pyrrolidine (PHP)	25	100
1-[1-(2-Thienyl)cyclohexyl]-piperidine (TCP)	50	50
1-(1-phenylcyclohexyl)-4-hydroxypiperidine	100	25
4-Phenyl-4-piperidinocyclohexanol	300	8.3
Ketamine	>500,000	<0.01

**1** Gently tear open the foil pouch to remove the stick. Pull the protective slide back to expose the sample pad.

**2** Immerse the sample pad into the specimen until the specimen reaches the DIP LINE. Hold the stick in the specimen for approximately 10 seconds.

**3** Remove the stick from the specimen. While holding the stick over the specimen cup, push the protective slide back over the sample pad.


Test complete—  
ready to interpret  
results

Test NOT complete—  
DO NOT interpret  
results

**4** Allow the test to proceed until a blue band forms in the “TEST VALID” window (see example on the LEFT). Test is NOT complete until a blue band forms in the “TEST VALID” window (see example on the RIGHT). Once the “TEST VALID” line appears, remove the cover tab and interpret results.

**5** A blue band in the result window indicates a NEGATIVE result.

**6** No color in the result window indicates a POSITIVE result. The window appears as a white plus (+) sign.

### Cross-Reactivity with Additional Drugs

The following compounds were tested at 100,000 ng/mL and found not to cross-react with the OnTrak TestStik for PCP assay.

Acetaminophen	Erythromycin
Acetylsalicylic acid	Estriol
Aminopyrine	17 $\alpha$ -Ethinylestradiol
Amitriptyline	Fenoprofen
Amobarbital	Furosemide
Amphetamine	Gentisic acid
Ampicillin	Glutethimide
Ascorbic acid	Guaiacol glycerol ether
Aspartame	Hydrochlorothiazide
Atropine	<i>p</i> -Hydroxyamphetamine
Benzocaine	3-Hydroxytyramine
Benzoylcegonine	Ibuprofen
Benzphetamine	Imipramine
Brompheniramine	Isoproterenol
Butabarbital	Ketamine
Caffeine	LSD
Calcium hypochlorite	Lidocaine
Chlordiazepoxide	MDA
Chloroquine	MDMA
Chlorpheniramine	Melanin
Chlorpromazine	Meperidine
Clemastine	Methadone
Cocaine	Methamphetamine
Codeine	Methapyrilene
Cyclizine	Methaqualone
Desipramine	Methylphenidate
Dextromethorphan	Methyprylon
Dextropropoxyphene	Morphine sulfate
Diazepam	Naloxone
Diphenhydramine	Naltrexone
Diphenylhydantoin	Naproxen
Dopamine	Niacinamide
Doxylamine	Nordiazepam
Ecgonine	Norethindrone
Ecgonine methyl ester	<i>l</i> -Norpseudoephedrine
<i>d,l</i> -Ephedrine	Nortriptyline
<i>d</i> -Ephedrine	Oxazepam
Epinephrine	Penicillin G

Pentobarbital	Quinine
$\beta$ -Phenethylamine	Secobarbital
Phenobarbital	Sulindac
Phenothiazine	Tetracycline
Phentermine	$\Delta^9$ -THC-9-carboxylic acid
Phenylbutazone	Tetrahydrozoline
Phenylephrine	Thioridazine
Phenylpropanolamine	Trifluoperazine
Phenyltoloxamine	Trihexyphenidyl
Procaine	Trimipramine
Procyclidine	Tripelennamine
Promethazine	Tyramine
Pseudoephedrine	Verapamil
Quinidine	Zomepirac

### References

1. Mandatory Guidelines for Federal Workplace Drug Testing Programs. *Federal Register* September 30, 1997. Department of Health and Human Services, Substance Abuse and Mental Health Services Administration.
2. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
3. Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc: 1984.
4. Baselt RC, Cravey RH. *Disposition of Toxic Drugs and Chemicals in Man*. 4th ed. Foster City, CA: Chemical Toxicology Institute; 1995:602-605.
5. Clouet DH, ed. Phencyclidine: An Update. *National Institute on Drug Abuse (NIDA) Research Monograph 64*. 1986.
6. Lundberg GD, Gupta RC, Montgomery SH. Phencyclidine: Patterns seen in street drug analysis. *Clin Tox*. 1976;9(4):503-511.
7. Lin DCK, Fentiman AF Jr, Foltz RL, Forney RD Jr, Sunshine I. Quantification of phencyclidine in body fluids by gas chromatography chemical ionization mass spectrometry and identification of two metabolites. *Biomed Mass Spec*. 1975;2:206-214.
8. Wong LK, Biemann K. Metabolites of phencyclidine. *Clin Tox*. 1976;9(4):583-591.
9. Data on file at Varian, Inc.

[www.varian-onsite.com](http://www.varian-onsite.com)

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#### Literatur

1. Mandatory Guidelines for Federal Workplace Drug Testing Programs. *Federal Register* 30. September 1997. Department of Health and Human Services, Substance Abuse and Mental Health Services Administration.
2. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
3. Blum K. *Handbook of Abusable Drugs*. 1<sup>st</sup> ed. New York, NY: Gardner Press Inc: 1984.
4. Baselt RC, Cravey RH. *Disposition of Toxic Drugs and Chemicals in Man*. 4<sup>th</sup> ed. Foster City, CA: Chemical Toxicology Institute; 1995:602-605.
5. Clouet DH, ed. Phencyclidine: An Update. *National Institute on Drug Abuse (NIDA) Research Monograph 64*. 1986.
6. Lundberg GD, Gupta RC, Montgomery SH. Phencyclidine: Patterns seen in street drug analysis. *Clin Tox*. 1976;9(4):503-511.
7. Lin DCK, Fentiman AF Jr, Foltz RL, Forney RD Jr, Sunshine I. Quantification of phencyclidine in body fluids by gas chromatography chemical ionization mass spectrometry and identification of two metabolites. *Biomed Mass Spec*. 1975;2:206-214.
8. Wong LK, Biemann K. Metabolites of phencyclidine. *Clin Tox*. 1976;9(4):583-591.
9. Dokumentation, Varian, Inc.

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### Bibliographie

1. Mandatory Guidelines for Federal Workplace Drug Testing Programs. *Federal Register* September 30, 1997. Department of Health and Human Services, Substance Abuse and Mental Health Services Administration.
2. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
3. Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc: 1984.
4. Baselt RC, Cravey RH. *Disposition of Toxic Drugs and Chemicals in Man*. 4th ed. Foster City, CA: Chemical Toxicology Institute; 1995:602-605.
5. Clouet DH, ed. Phencyclidine: An Update. *National Institute on Drug Abuse (NIDA) Research Monograph 64*. 1986.
6. Lundberg GD, Gupta RC, Montgomery SH. Phencyclidine: Patterns seen in street drug analysis. *Clin Tox*. 1976;9(4):503-511.
7. Lin DCK, Fentiman AF Jr, Foltz RL, Forney RD Jr, Sunshine I. Quantification of phencyclidine in body fluids by gas chromatography chemical ionization mass spectrometry and identification of two metabolites. *Biomed Mass Spec*. 1975;2:206-214.
8. Wong LK, Biemann K. Metabolites of phencyclidine. *Clin Tox*. 1976;9(4):583-591.
9. Documentation de Varian, Inc.

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1. Mandatory Guidelines for Federal Workplace Drug Testing Programs. *Federal Register* 30 septiembre 1997. Department of Health and Human Services, Substance Abuse and Mental Health Services Administration.
2. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
3. Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc: 1984.
4. Baselt RC, Cravey RH. *Disposition of Toxic Drugs and Chemicals in Man*. 4th ed. Foster City, CA: Chemical Toxicology Institute; 1995:602-605.
5. Clouet DH, ed. Phencyclidine: An Update. *National Institute on Drug Abuse (NIDA) Research Monograph 64*. 1986.
6. Lundberg GD, Gupta RC, Montgomery SH. Phencyclidine: Patterns seen in street drug analysis. *Clin Tox*. 1976;9(4):503-511.
7. Lin DCK, Fentiman AF Jr, Foltz RL, Forney RD Jr, Sunshine I. Quantification of phencyclidine in body fluids by gas chromatography chemical ionization mass spectrometry and identification of two metabolites. *Biomed Mass Spec*. 1975;2:206-214.
8. Wong LK, Biemann K. Metabolites of phencyclidine. *Clin Tox*. 1976;9(4):583-591.
9. Datos de archivo, Varian, Inc.

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## Riferimenti

1. Mandatory Guidelines for Federal Workplace Drug Testing Programs. *Federal Register* September 30, 1997. Department of Health and Human Services, Substance Abuse and Mental Health Services Administration.
2. Hawks RL, Chiang CN (a cura di) Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
3. Blum K. *Handbook of Abusable Drugs*. prima ed. New York, NY: Gardner Press Inc: 1984.
4. Baselt RC, Cravey RH. *Disposition of Toxic Drugs and Chemicals in Man*. quarta ed. Foster City, CA: Chemical Toxicology Institute; 1995:602-605.
5. Clouet DH (a cura di) Phencyclidine: An Update. *National Institute on Drug Abuse (NIDA) Research Monograph 64*. 1986.
6. Lundberg GD, Gupta RC, Montgomery SH. Phencyclidine: Patterns seen in street drug analysis. *Clin Tox*. 1976;9(4):503-511.
7. Lin DCK, Fentiman AF Jr, Foltz RL, Forney RD Jr, Sunshine I. Quantification of phencyclidine in body fluids by gas chromatography chemical ionization mass spectrometry and identification of two metabolites. *Biomed Mass Spec*. 1975;2:206-214.
8. Wong LK, Biemann K. Metabolites of phencyclidine. *Clin Tox*. 1976;9(4):583-591.
9. Dati presenti nel file, Varian, Inc.

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### Bibliografia

1. Mandatory Guidelines for Federal Workplace Drug Testing Programs. *Federal Register* September 30, 1997. Department of Health and Human Services, Substance Abuse and Mental Health Services Administration.
2. Hawks RL, Chiang CN, eds. Urine testing for drugs of abuse. *National Institute on Drug Abuse (NIDA) Research Monograph 73*. 1986.
3. Blum K. *Handbook of Abusable Drugs*. 1st ed. New York, NY: Gardner Press Inc: 1984.
4. Baselt RC, Cravey RH. *Disposition of Toxic Drugs and Chemicals in Man*. 4th ed. Foster City, CA: Chemical Toxicology Institute; 1995:602-605.
5. Clouet DH, ed. Phencyclidine: An Update. *National Institute on Drug Abuse (NIDA) Research Monograph 64*. 1986.
6. Lundberg GD, Gupta RC, Montgomery SH. Phencyclidine: Patterns seen in street drug analysis. *Clin Tox*. 1976;9(4):503-511.
7. Lin DCK, Fentiman AF Jr, Foltz RL, Forney RD Jr, Sunshine I. Quantification of phencyclidine in body fluids by gas chromatography chemical ionization mass spectrometry and identification of two metabolites. *Biomed Mass Spec*. 1975;2:206-214.
8. Wong LK, Biemann K. Metabolites of phencyclidine. *Clin Tox*. 1976;9(4):583-591.
9. Ficheiros da Varian, Inc.

www.varian-onsite.com

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Dentro dos E.U.A.: 1-800-737-9667

Fora dos E.U.A.: 00 1 (949) 770-9381, extensão 2

DAT.TechSupport@varianinc.com

 Varian, Inc.  
Lake Forest, CA 92630 USA

 Varian B.V.  
Herculesweg 8, 4338 PL Middelburg  
 Holanda

Número da patente: 5,770,458

16153 1

© 2003, Varian, Inc.

**VARIAN**